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1. Table of contents

1. Table of contents	60
2. Introduction	62
2.1 About these instructions for use.....	62
2.2 Symbols used in the text.....	62
2.3 Intended use.....	63
2.4 General safety notes	64
3. Product description and configuration	65
3.1 Product description	65
3.2 Configuration.....	66
4. Start-up.....	67
4.1 Delivery	67
4.2 Cable version connection to the nurse call system.....	67
4.3 Wireless version connection to the nurse call system.....	68
4.4 Multi-call function (only available with wireless version).....	69
4.5 Operating aids/Information videos	70
5. Operation.....	71
5.1 Setting up the low nursing bed	71
5.2 CMC-Box (ConnectivityMultiCall-Box)	72
5.3 Activate the battery box.....	72
5.4 Battery operation.....	73
5.5 Moving and braking (sentida 7-i cozy)	73
5.6 Moving and braking (sentida 7-i mobile)	74
5.7 Calf section mechanically adjustable	76
5.8 Emergency lowering of backrest	76
5.9 Adjust SafeFree® side guards	77
5.10 Mobi-Lift®	79
5.11 Height adjustment.....	80
5.12 Handset.....	81
6. CareBoard®.....	83
6.1 Positioning.....	85
6.2 Locking functions	86
6.3 Safety sensors	87
6.4 Bed exit assistant.....	89
6.5 Weighing system.....	90
6.6 Weight monitoring	91
6.7 Weight history	92
6.8 Nursing routine documentation.....	93

6.9 Nurse settings	94
6.10 Underbed light (UBL).....	95
6.11 Digital user manual.....	95
6.12 Service reminder (pop up).....	96
7. Accessories.....	96
7.1 Mattress	96
7.2 Lifting pole	97
7.3 Infusion stand	98
7.4 Extensions of patient surface	99
7.5 Foldable linen holder.....	100
7.6 Integrated linen holder	101
7.7 Accessories holder.....	101
7.8 Telescopic wall spacers.....	102
7.9 Side guard padding, side guard net	102
7.10 Extension for SafeFree® side guard (accessory).....	103
8. Cleaning and disinfection.....	104
8.1 General procedure for cleaning/disinfection.....	104
8.2 Cleaning.....	105
8.3 Disinfection	105
8.4 Special information regarding real wood parts.....	105
9. Maintenance and service.....	106
9.1 Spare parts	107
9.2 Troubleshooting.....	108
9.3 Storage.....	110
10. Disposal.....	111
11. Warranty	111
12. Technical data	112
13. Classification.....	114
13.1 Applied standards	114
14. Product identification	115
14.1 Symbols used	115
14.2 Nameplate.....	116
14.3 Scale nameplate	117
15. Declaration of Conformity.....	351

2. Introduction

2.1 About these instructions for use

In this section, you will find information about the structure of these instructions for use and a description of the signs and symbols used. This document contains instructions for use of the digital nursing bed sentida 7-i.

These instructions for use may contain inaccuracies and typographical errors. The information provided in these instructions for use are updated at regular intervals, and changes resulting from product maintenance are implemented in future editions. Changes or improvements are possible at any time without previous notice. Contact our customer service in case you have any questions.

The instructions for use have to be read and observed by every person operating the digital nursing bed.

Apart from the instructions for use and the accident prevention regulations valid for the respective country and area of use, the commonly accepted regulations for safe and professional work must also be adhered to.

Along with these instructions for use, you have the possibility to use the brief, digital user manual, which is integrated in the CareBoard®. There, you can find a short overview of the nursing home bed's components, mechanical adjustment options and detailed notes about the use of the CareBoard®.

2.2 Symbols used in the text

In these instructions for use, the following terms and symbols are used for important notes:

Danger!

Safety notes on people being in danger are marked with this symbol. This symbol is used when there is imminent danger of death or severe injuries.

Attention!

This symbol is used when there might be the danger of unsafe situations or minor injuries.

Caution!

This symbol precedes warnings when there is the danger of damage to the equipment or other things.

This symbol precedes additional helpful advice.

- A hyphen preceding the text means: This is part of an itemization.
- A dot preceding the text means: You have to adhere to this instruction.
Indented text means that this is the result of your actions.

2.3 Intended use

The digital nursing bed sentida 7-i is a Class Im (non-sterile with measuring function) medical device. It is intended exclusively for human medical treatment and should only be used in the application environments 3 and 5 in accordance with the standard IEC 60601-2-52:2010 + A1:2015: Application environment 3 includes longterm care in nursing homes and homes for the elderly, in geriatric and rehabilitation facilities. These institutions require medical supervision and, if necessary, monitoring. A ME DEVICE used for medical purposes can be made available to maintain or improve the resident's condition.

Application environment 5 includes the ambulant care area in a hospital or other medical institution with medical supervision at hand. A ME DEVICE is made available for persons with illness, injuries or disabilities for treatment, diagnosis or supervision. Any other application must be agreed beforehand in writing with Wissner-Bosserhoff GmbH. The product must be used as nursing equipment and is subject to the regulations of the relevant professional organization.

This nursing bed is a medical product in the sense of the regulations and standards named in the following chapter „product safety / excerpt of applied standards“. Accordingly, this product may only be used under medical supervision. The decisive factor in whether nursing is medically supervised or not depends on whether nursing takes place under the instructions of medical personnel or not.

It is designed for a safe working load of 250 kg (lifting force), with a maximum patient or resident weight of 215 kg and a mattress and accessory weight of 35 kg. If other, heavier or further loose or attached accessories are used, the maximum patient weight is reduced accordingly.

The bed may only be operated by persons who can guarantee correct handling because of their training or knowledge and practical experience. The users of the digital nursing bed must have received instructions for correct use and must have familiarized themselves with the equipment with the help of these instructions. The digital nursing bed may only be used in complete accordance with the instructions for use.

We assume no liability for potential product damage or personal injuries caused by unapproved accessories or mutual cancellation of the intended use.

The specified use is the intended use. Its sole source for the operator or user is the labeling and the instructions for use. The digital nursing bed sentida 7-i may only be used in full accordance with the instructions for use.

The bed can be moved and therefore the easy positioning within the room ensures easy cleaning of the digital nursing bed and the space underneath. It is, however, not suitable for the transport of residents and patients (only sentida 7-i | cozy).

Please note that the system is not designed to be used as an emergency call or a lifesaving call for help. It is rather a product that helps to facilitate the daily work in a nursing home.

2.4 General safety notes

The digital nursing bed sentida 7-i was designed in accordance with state-of-the-art technology and the approved safety regulations. Nevertheless, its use can endanger the life or health of the user or third parties and have adverse effects on the digital nursing bed sentida 7-i and other material assets.

Only use the digital nursing bed sentida 7-i when it is in perfect working order. Only use it in compliance with its intended use, safely and aware of dangers and in accordance with these instructions for use! Malfunctions that may affect safety must be eliminated immediately in particular!

Always keep these instructions for use ready to hand at the digital nursing bed's operating place. In addition to the instructions for use, please observe the generally applicable standards and commonly accepted regulations on accident prevention and environmental protection!

Without permission of the manufacturer, do not make any changes or conversions to the digital nursing bed sentida 7-i. Spare parts have to comply with the manufacturer's requirements. For original spare parts, this is always guaranteed.

Adhere to the prescribed checks and tests! Make sure that operating and auxiliary materials and replacement parts are recycled safely and in an environment-friendly manner.

Before operating the bed, ensure the functional safety and the proper condition of the bed. Always leave the castors braked; this avoids the risk of falls for the patient/resident when getting in or out of the bed.

Adjust the height of the patient surface to the size of the patient to avoid the risk of falls. If incontinent patients occupy the bed, incontinence protection must be used for the mattress.


3. Product description and configuration

In the following, the digital nursing bed sentida 7-i is also referred to as nursing bed or bed.

The terms left and right have to be seen from the perspective of the person lying in bed.

3.1 Product description

- | | |
|--------------------------------------|--|
| 01 Fußteil | 08 Undercarriage cladding |
| 02 Mobi-Lift® | 09 Brake lever |
| 03 Calf section | 10 Thigh section |
| 04 Socket for humidity sensor | 11 Side panel |
| 05 Backrest | 12 Comfort resident handset SafeControl |
| 06 CareBoard® | |
| 07 Head section | |

3.2 Configuration

Depending on the configuration chosen, the sentida 7-i can have different functions and features:

Standard version:

Weight monitoring

- integrated weighing system
- weight measurement at the push of a button

BMI check

- identifying malnutrition
- calibrated scale

Bed exit manager

- early warning at edge of bed
- timing of nurse call can be set as required

Connectivity

- integrated interface standards
- transmission in real time

Multicall (only with wireless version)

- calls are classified
- prioritization possible

System security

- closed system
- secure access management

Additional options:

Nursing documentation

- paperless
- data entry directly by the bed

Moisture kit

- detects moisture early
- choice of 2 sensor mats

Safety sensor

- early warning at edge of bed
- timing of nurse call can be set as required

Mobilization kit

- light: plug & play handset
- full: light feature + Mobi-Lift®

4. Start-up

4.1 Delivery

The nursing home bed is generally completely assembled when delivered, or is assembled by specialists on site.


- Please check whether all components have been delivered by consulting the delivery papers.
- Write down potential defects or damage on the delivery note.
- Inform your service partner immediately about potential transport damage or defects. The address and phone number are noted on the last page of this manual.

➔ Upon delivery, the time of the CareBoard® needs to be set depending on the time zone (refer to section 6.11 for details).

4.2 Cable version connection to the nurse call system

If the CareBoard® is attached at the head-end (e.g. wall mounting), the following cable exits are in the center below the head end (see picture):

- a** Network cable connection (network cable not provided)
- b** NurseCall cable
- c** Cable CareBoard®
- d** Adapter cable NurseCall (nurse call system specific) with or without ON/OFF switch
- e** Power cable

In case of a wall mounting, first mount the enclosed CareBoard® holder to an easily accessible place at a user-friendly height.

Attention!

If the wall mounting including the CareBoard® comes loose, injuries could be caused. To avoid this, we recommend a safe wall anchoring next to the nursing bed. Mounting material (dowels and screws) for wall mounting is not enclosed, as it depends on the consistency of the wall. Please use mounting material that is suitable for the respective consistency of the wall. In the case of related questions, please contact the local specialist shop.

If the CareBoard® is attached to the foot section of the bed with a bracket (optional), this bracket needs to be mounted about 27 cm from the lower edge of the foot section. If the bed has an „A-k“ foot section, this spacing is reduced to about 15 cm. Pass the CareBoard® cable out centrally under the foot section.

If the CareBoard® is attached to the bed's foot end with a holder (optional), the CareBoard® cable exit is in the center below the foot section.

The adapter cable (pos. d) must be laid out individually for the call systems installed in the building and must be ordered with Wissner-Bosserhoff GmbH. Connect the bed to the nurse call system by attaching the nurse call cable (pos. b) to the adapter and the other end to the call system.


To allow you to continue using the hand-held switch on your nurse call system and also to connect sentida 7-i, two sockets are necessary. If your nurse call system does not have a socket, a Y adapter is required. You can obtain one from your call system manufacturer. Please also ask the manufacturer about any risks and other information that should be noted.

If there is a Y cable for the corresponding call system, both connectors on the Y cable may need to be occupied (e.g. by a handheld switch). As final step, connect the power cable to the power supply.

Wissner-Bosserhoff is not responsible for the incorrect use of Y adapters or any possible damage to the call system that may result, and accepts no liability for such damage.

4.3 Wireless version connection to the nurse call system

The wireless version connection for the most part corresponds to the cable version connection. The only difference is that no cable is connected to the nurse call connector and an adapter cable is connected to the wireless receiver instead (see following illustration). The only cables coming from the bed are the CareBoard® and power cable. With regard to the previous illustration, this means that position c is omitted.




→ For optimal signal transmission, the control unit should not be placed near any large metal surfaces. If possible, it should be attached to wooden or plastic elements.

4.4 Multicall function (only available with wireless version)

The multicall function makes it possible to distinguish between different safety calls activated at the same time and display these on a compatible end device for communication purposes. As a result, nursing staff are able to prioritize calls and therefore decide which residents need help most urgently. This is only possible if the nurse call system in question is able to read and display individual calls.

5 different safety calls

Calls can be displayed individually on PC, DECT, or smartphone!

4.5 Operating aids/Information videos

The following QR codes can be scanned to call up information videos on the various features of the sentida 7-i:

BMI

Nursing
documentation

Scale

Multicall

System security

5. Operation

5.1 Setting up the digital nursing bed

For setting up the nursing bed sentida 7-i, rooms with even hard floors are suitable.

Floors that are, for example, too soft, unsealed or with imperfections are not suitable. Soft wooden floors, porous and soft stone floors, carpets with foam backing, soft linoleum or similar floor covering are mostly unsuitable.

Please make sure that the bed is positioned freely in the room. When setting up the bed in the intended room, there has to be a safety space around the bed of 20 cm to the wall, window sills or other furniture. In each position of the bed, also when adjusting the height, there has to be a distance of 2,5 cm to the wall, other objects and installations.


- When choosing a location, check for mechanical and electrical protection measures such as wall-mounted bumpers, corner guards, RCB switches and the like.
- Pay attention to the whereabouts of suitable wall sockets for the power supply connection for the bed. Positions that could lead to collisions or stretching when the bed is adjusted, or in which the power supply cable stretches too far across the floor or even under the bed, are unsuitable.
- Extension cords and multiple adapters which lie unsecured on the floor must not be used.

➡ The nursing home bed has to be positioned in the room in a way that it can be easily disconnected from the power supply.

➡ The nursing home bed is usually not equipped with an equipotential connection. A combination with other electro-medical equipment is therefore not permitted without further safety measures. An equipotential connection can be optionally provided, if required. For further information, get in touch with the technical customer service.

5.2 CMC-Box (ConnectivityMultiCall-Box)

→ The CMC box can be found at the head end on the left under the side beam and only needs to be used if a reset is required.

a Reset-button

b Network cable

c Nursecall cable


d Connection cable Powerbox

5.3 Activate the battery box

Attention!

When removing the insulating foil from the battery box, please be careful as to not sustain cutting injuries. Wearing protective gloves is recommended.

Battery box
Insulating foil (red)

Removing the insulating foil

- Remove the red insulating foil from the battery box to clear the charging pins.
- Check whether the insulating foil is complete and intact as shown below.

- Insulating foil (red)
- In the event of damaged insulating foil, please contact the Wissner Bosserhoff customer service (see last page).

5.4 Battery operation

The bed is equipped with a battery, which means it can be adjusted independently from the supply network.

The number of adjustments depends on the charge status of the battery. Complete discharging of the battery significantly reduces its service life. To prevent the batteries from excessive discharge, the bed should remain connected to the power supply. The battery charges automatically via an integrated charging board

- The battery is maintenance-free
- It is integrated in the power box
- It is a lead gel battery which can be operated irrespective of its location.
- Nominal data: 24V, 1.2Ah
- Depending on the usage, the service life of the battery is approx. 5 years
- Disposal is only permitted at special collection points.

Attention!

Harmful chemicals can leak from opened batteries and cause injuries. Never open the battery box.

Exchanging the batteries may only be performed by trained personnel. The service life of the battery strongly depends on its use. Complete discharging of the battery significantly reduces its service life.

5.5 Moving and braking (sentida 7-i | cozy)


The bed can be moved and therefore the easy positioning within the room ensures easy cleaning of the nursing bed and the space underneath. However, the nursing bed is not suitable for the transport of residents and patients. The nursing bed may only be moved in its lowest position and without a patient. Moving the nursing bed to other rooms with a resident or patient still in it does not comply with the intended use.

There is a covered axle with lockable castors beneath the head and foot end. On the inner side of this undercarriage cover, there is a braking lever located at the left and right side of the bed. With this lever, you can release or lock the castors for each axle.

Caution!

The bed can only be moved if both brakes are released. If braked castors are moved across the floor, the flooring could be damaged.

1. Pull the braking lever up to release the brakes. Pictured in released position.
2. Push the braking lever down to lock the nursing bed. Pictured in locked position.


5.6 Moving and braking (sentida 7-i mobile)

To activate the brake bar, move the bed into the bed egress position or higher. The brake bar can now be activated with the foot in three positions.

Position 01: Brake

Position 02: Can be moved

Position 03: Direction locking
(The tilting movement of the left castor at the head end is blocked, thereby making it easier to maneuver the bed.)

When moving the bed, please adhere to the following steps:

- Release the brakes.
- At its head or foot end, move the bed to the desired position.
- Then, activate the brakes again by pushing the brake bar .

**Caution!**

The side panels' anchoring could come loose under high load. Only move the bed at its head or foot end.

**Attention!**

Injuries could be caused if feet or other body parts get caught in the moving castors or the undercarriage cladding when moving the bed. When moving the bed, make sure that your foot does not get caught in the moving castors or the undercarriage cladding.

**Caution!**

Power supply connection cables could break when they are run over. Never run over power supply cables or other cables.


The nursing bed sentida 7-i is equipped with a spiral cable. A plastic hook is attached to the spiral cable.

- If you want to move the nursing bed sentida 7-i, hang the power supply connection cable's hook to the upper edge of the head end. Thus, you can avoid that the power supply connection cable is run over.
- Regularly check the power supply connection cable, the correct insertion of the power plug, the power plug cover and the housing openings on the bed for signs of damage, in order to rule out hazards.

5.7 Calf section mechanically adjustable


The height of the calf section can be adjusted by using ratchets.

- If you want to lift the calf section manually, lift it on the frame. It clicks into the desired position.
- If you want to lower the calf section manually, lift it on the frame as high as possible, then lower it.

5.8 Emergency lowering of backrest

The backrest can be quickly lowered manually, without handset and CareBoard®. For this purpose, a release handle with a Bowden cable is situated at the head end, below the backrest. By using the release handle, uncouple the backrest motor with the Bowden cable. Afterwards, the released backrest can be lowered manually.

- Lift the mattress above the backrest.

The release handle is now accessible.

- Hold the handle on the backrest frame with one hand. Squeeze the release handle with the other hand while pressing the backrest against the weight of both, the frame and the resident

With the help of the Bowden cable, the backrest motor is decoupled and the backrest is pushed down by its weight.

- Slowly lower the backrest.

➔ The manual emergency lowering of backrest is intended for emergency use only. Please use only when necessary.

5.9 Adjust SafeFree® side guards

When using the side guards, particular care and caution should be taken: Side guards are first and foremost intended to protect the patient. On no account may the side guards be used as a means to restrain the patient.

➔ As soon as the mobility of the patient or resident is restricted, a court order is required to use the side guards. Divided side guards or side guards that can be individually lowered, in which only the head and section is assembled, and which can be used as a help in getting in and out of bed, are only a permissible exception for mobile residents or patients.

➔ When using side guards, their suitability should be checked, taking into account the individual requirements of the respective patient or resident. In particular, the clearances between the panels and supports in relation to the body shape of the patient or resident must be taken into account, so that entrapment or slipping through is ruled out. If necessary, additional side guard cushioning must be used which reduces the size of the gap or closes it. This decision is exclusively the responsibility of the medical supervision in charge.

Attention!

If necessary, use a side guard extension to increase the height when using raised mattress systems. Also use the side guard cushioning available for this (accessories) as an important additional safety and protection feature.

Attention!

When lowering or raising the side guards, fingers, hands, or other body parts of the resident could be squashed. Raise or lower the side guards with the highest care. Only hold the side guards at the positions indicated in the instructions for use.

Attention!

Restless and disoriented people can become trapped in the free spaces of the side guard. Protect these by using side guard padding.

**Attention!**

Side guards not locked correctly can slide down and cause the resident to fall. By shaking the side guard, make sure that it has clicked into the right position. The correct engagement can also be visually inspected: The yellow locking buttons in the locking recess of the telescope mechanism must be recognizable on both sides.

**Caution!**

The SafeFree® side guard could get overloaded and damaged by heavy loads, leaning or pulling. Such overloading has to be avoided.

At each side of the bed, the SafeFree® side guard is divided into two parts. There is a release button on the bottom side of both ends of the handrail, with which the lock can be unlocked. It can be extended in 4 steps up to the height of 400 mm.

→ Move the side guard with both hands to ensure smooth lowering/raising in both guide rails.

Release buttons

1. level: Barrier-free

Side guard is in the lowered position.

**2. level: Mobilization**

3/4 protection at a height of 242 mm.

**3. level: Safety**

Complete protection at a height of 340 mm.

**4. level: Increased safety**

More safety at a height of 400 mm for thicker mattresses.

In case you want to raise the SafeFree® side guard, follow these instructions:


- Grab the SafeFree® side guard at its handrail; lift it, until the SafeFree® side guard audibly clicks into the next position.
- To further lift the SafeFree® side guard, press both release buttons simultaneously and pull the side guard upwards at its handrail until the SafeFree® side guard audibly clicks into the next position.

To lower the SafeFree® side guard, proceed in reverse order.

- Press both release buttons **simultaneously** and lower the SafeFree® side guard at its handrail to the desired position.

5.10 Mobi-Lift®

The low nursing bed sentida 7-i is equipped with a Mobi-Lift®. The resident uses Mobi-Lift® for bed egress and ingress, as the bed can be lifted or lowered.

Mobi-Lift® in horizontal (left) and in fixed position (right)

The Mobi-Lift® is attached to the left and right of the calf section's frame. Here, it can be pulled up and fixed, or positioned horizontally.

- If you want to pull up Mobi-Lift®, lift it from its horizontal to a vertical position and fix it there by letting it glide into its holder.
- In order to put aside the Mobi-Lift®, vertically pull it out of its holder and place it in a horizontal position.

The Mobi-Lift®'s handle is equipped with two push buttons with which the bed can be lifted or lowered.

When the SafeFree® side guard is lowered, a mobile resident can set the height of the bed by using the two push buttons. As soon as a safe and comfortable height to leave the bed is reached, the resident can leave or enter the bed.

5.11 Height adjustment

Danger!

Before every height adjustment of the bed, you must make sure that there are no persons, animals or objects in the danger area. Always leave the handset in the locked position so that no unauthorized persons can manipulate the height adjustment of the bed.

Attention!

When lowering the bed or tilting it to the seating position, make sure not to get your feet in the danger area (at the bed ends or under the bed).

Attention!

Comfort seating position / legs up position: Particularly in conjunction with a bed extension, there is a danger of entrapment between the extension and the floor – and of the foot section or head section coming into contact with the floor. Move the patient surface to an appropriate height before every adjustment of the bed. The same applies when a linen holder is used!

Attention!

Due to its design as a low nursing bed, there is a clearance of less than 150 mm. Prior to every height adjustment of the bed, you must make sure that there are no persons, objects or anything else in the entrapment area beneath the bed! Always leave the handset in the locked position so that no unauthorized persons may adjust the height of the bed.

Bed egress position

When adjusting the height of the bed upward or downward, the bed stops at the comfort bed egress position. To continue adjusting the bed, the function key on the handset or CareBoard® needs to be released and then pressed again. Using movable lifting equipment (e. g. mobile patient lift)


The low nursing bed's design requires a safety clearance of less than 150 mm. Make sure before each use of movable lifting equipment that the free space beneath the low nursing bed is sufficient for the movable lifting equipment. In any other case, adjust the height before using movable lifting equipment.

5.12 Handset

The comfort resident handset is also called handset in this manual.

The handset is included in the delivery of the low nursing bed sentida 7-i. It is held by a holder for the comfort resident handset. You can attach a bracket for the comfort resident handset to the SafeFree® side guard on the left or right at the head end (accessory).


- If the handset is to be switched from the left SafeFree® side guard to the right one, pull the multiple pin plug from the left connection and seal the socket with the attached stopper.
- Connect the plug with the right connection after the stopper has been removed from the socket.
- Guide the connection cable together with the handset of the right connection past the undercarriage cover underneath the right side panel and up the outer side of the side panel to the handset holder.
- The buttons of the handset are illuminated when the bed is connected to the power supply.

➔ Individual functions or all functions of the handset can be locked. This can be recognized by the illuminated LEDs of the handset.

**Caution!**

The comfort resident handset is sensitive to mechanical forces and liquids. Avoid that the handset falls to the ground or gets overloaded. Avoid the penetration of liquids into the handset. Clean the handset with cleaning cloths for computer keyboards. Regularly check the handset and the cables for signs of damage. Defective handsets have to be replaced immediately.

→ The use of the handset by the resident/patient himself lies exclusively at the discretion and responsibility of the medical supervisor and nursing staff in charge. Corresponding measures should be noted in the shift report or in the nursing documentation to guarantee a correct shift handover and traceability of usage rights.

Before you adjust a position of the low nursing bed sentida 7-i with the handset, you have to activate the handset with the GO safety button. It is the green button with the button symbol.

- Activate the handset with the GO safety button

The handset is now activated. It automatically switches off after 3 minutes without having to press the button.

As long as the membrane keyboard is illuminated, the handset is activated and you can adjust all released electrical functions.

In the upper gray area, you can adjust the leg rest on the left and the backrest on the right. In the orange-colored field below, you can adjust the height on the left and the auto contour on the right.


- Press the corresponding button in the gray or orange-colored field to adjust a position.
- Deactivate the handset on the CareBoard® after you have adjusted the required position. This way, you avoid unintentional misuse.

There is a LED light integrated into the handset.

- If you want to switch on the LED light on the handset, press the corresponding button in the lower area of the handset.

6. CareBoard®

The CareBoard® provides access to all technical adjustment options of the low nursing bed.

**Danger!**

Liquids on the CareBoard[®] display can lead to injuries caused by movements of the low nursing bed. After use, always press the emergency stop button. Immediately wipe liquids off the display and also immediately dry it after cleaning.

**Caution!**

When mounting the CareBoard[®] to the wall, please make sure that it cannot collide with the bed when it is moved or adjusted.

If the CareBoard[®] is separated from the bed, the handset also ceases to function.

If 10 minutes pass without any action via the CareBoard[®], the user is automatically logged out and the equipment switches to standby mode. This can also be performed manually by pressing the emergency stop button.

There is a bank of buttons to the left of the display. At the top, the bank of buttons shows the GO safety button which activates the CareBoard[®]; that is why this button is called GO button. There are 3 further buttons below it, which offer access to important acute functions. This direct access is also called acute care pack.

Difference between display off and log out: When the CareBoard[®] is activated, pressing the GO safety button once switches off the display. The user is only logged out if the emergency stop button is pressed.

You can attach a bracket for the comfort resident handset to the SafeFree[®] side guard on the left or right at the head end.

This acute care pack contains the following functions:

- Automatic electrical CPR

CPR (cardiopulmonary resuscitation), this corresponds to the comfort bed egress position with horizontal position of knee and backrest. The CPR automatic is used during reanimation by cardiac massage.

- Automatic electric Trendelenburg

Here, you can set the Trendelenburg position. The patient surface is even and the patient's head is placed below the patient's center of gravity on a surface which is inclined at an angle of at least 12°.

- Emergency stop button

By pressing the emergency stop button, all electric motor adjustments are stopped immediately and the CareBoard[®] is switched off.

The symbols below the bank of buttons indicate if the low nursing bed is connected to the power supply system and if the battery is loading. The following table shows the meanings of the battery symbol.

The battery symbol flashes.	The battery is loading.
The battery symbol lights up permanently.	The battery is almost empty; it has to be charged; only a few adjustments can still be made.
The battery symbol remains dark.	The battery is fully charged.

- Activate the CareBoard® with the Safety-GO button.


The CareBoard® is now activated; however, it automatically switches off after 3 minutes without having to press the button.

The display is activated as long as it lights up, and you can adjust the low nursing bed sentida 7-i by electric motor. The adjustment possibilities are described on the following pages.

- Deactivate the CareBoard® by pressing the emergency stop button after you set the requested position. This way, you avoid unintentional misuse.

6.1 Positioning

In this menu, you can adjust the position of the bed by pressing the arrow buttons. The arrow buttons are normally displayed in black; however, the color changes to white when the indicator is being pushed. Only when this function is locked, the lock symbol is displayed.

01 Adjusting the leg rest

02 Auto contour (simultaneously adjusting the leg and backrest)

03 Adjusting the backrest

04 Height adjustment

05 Head down position

06 Head up position


07 Comfort seating position

08 Prevention position (also Semi-Fowler position for bedsore prevention)

6.2 Locking functions

On the CareBoard®, locked functions are displayed with a yellow highlighted lock symbol; the locked functions cannot be carried out with the handset anymore. In the position menu, an additional lock symbol is displayed above the locked functions. If someone still tries to use a locked function, the CareBoard® temporarily shows a lock symbol and a signal sounds.

Besides single functions, you can also lock all functions simultaneously for operation via handset. No functions can be used via the handset, until they are unlocked on the CareBoard®.


- 01** Locking the handset
- 02** Locking the backrest
- 03** Locking the leg rest
- 04** Locking the height adjustment

➔ If single functions are locked, combined functions like comfort seating or auto contour are carried out. The locked power units are not operated. Normally, this does not lead to the desired position. However, the functions of the acute care package are carried out without taking possible lockings into consideration.

➔ If the bed is not equipped with a handset, the latter is not displayed in the overview on the screen.

6.3 Safety sensors

A red highlighted, flashing speaker symbol indicates a current light alarm. A green speaker symbol indicates the proper condition of the monitored function. A speaker symbol that is not highlighted indicates that the function is not monitored.

01 Bed exit assistant: Submenu for switching to calls triggered at edge of bed

02 Humidity assistant: Submenu for sensitivity adjustment

03 Side guard assistant: Submenu with selectable modes for call triggering

04 Bed height assistant: Triggers a call when the bed height is changed

05 Nurse present button: In this mode, the nurse call system does not send any calls. It switches off automatically after around 15 minutes. Live calls are also acknowledged via this.

06 Brake assistant: Triggers a call, when a brake is released

07 Quit nurse call: All active (red marked) calls are reset

➡ The functions shown on the screen depend on the selected bed configuration.

➡ To access additional settings for individual functions (submenus), click the gear symbol (see next page for description).


A red highlighted, flashing speaker symbol in the CareBoard® display indicates a current light alarm. If the bed is connected to a call system, the light alarm is transmitted optically and acoustically depending on the call system. The light alarm can then be turned off at the call system.

- Now you have to find the cause of the light alarm on the low nursing bed sentida 7-i.
Example: A brake was released. Lock the brake once again.
- Then, turn off the light alarm in the alarm menu. To do so, press the crossed out speaker symbol.

The red highlighted flashing speaker symbol is now turning green. This indicates the proper condition of the monitored function.


➔ If the light alarm cannot be transmitted to the call system despite a properly working CareBoard®, the intermediate switch on the adapter cable (if used) might be switched „OFF“. Switch „ON“ the intermediate switch. If still no light alarm can be transmitted, a technical defect on the intermediate switch, the NurseCall box or the connection cable might be the case. In this event, please contact the wi-bo customer service.

➔ If a humidity sensor mat was not connected, the pictogram is crossed out in red.

Humidity assistant submenu

In this submenu, which you can open by pressing the gear symbol next to the humidity assistant, you can adjust the sensitivity of the humidity sensors from 1 to 8 (8 = maximum sensitivity).

Side guard assistant submenu


The submenu for the side guard assistant can also be opened by pressing the corresponding gear symbol. Here, you can select between „Safe“ and „Free“ mode. To switch modes, you have to enter a password in the field below.

„Safe“ mode: A call is triggered if any one of four raised side guards is lowered.

„Free“ mode: A call is triggered if the fourth side guard is raised when three side guards were already up.

6.4 Bed exit assistant

Delay times can be set for the light alarm when someone leaves the bed. During these tolerance times, the light alarm is suppressed. This enables mobile residents to go to the toilet without triggering the light alarm. The menu for setting the tolerance times can be opened by touching the wrench symbol in the previous menu.

- 01** You can set a delay time for this, during which the light alarm is suppressed while the resident is away from the bed. To cover visits to the restroom, for example, a tolerance time of 0 – 30 minutes is typical. If the resident does not return to the bed within the set time, a light alarm is triggered.
- 02** This selection determines whether or not a call is triggered if the resident has left the bed (bed exit) or is just getting up (early-warning alarm).
- 03** If the early-warning alarm option is selected, you can specify the call sensitivity in this sub-menu. Please also consider the resident's weight when setting this value.
- 04** This setting determines whether the bed exit assistant is active permanently or only in a definable period.
- 05** Set the period of time when the bed exit assistant is inactive.

➔ As the bed exit function is based on the weighing system, the mentioned settings can only be made while the resident is in bed. Otherwise, the system cannot transmit reliable signals.


6.5 Weighing system

In order to guarantee exact measuring results, place the bed on an even, horizontal ground. For checking purposes, please use the circular level which is attached below the backrest.

To measure, bring the bed back to a horizontal position. We recommend the position using the automatic CPR function („CPR“ button of the acute care pack on the CareBoard®).

In doing so, please make sure that the bed is positioned freely in the room and that the complete patient surface is not loaded or unloaded by external influences.

In the event of consecutive measurings, always bring the bed to an identical position in order to avoid measuring inaccuracies.


The circular level is attached in the center of the head-end cross bar

- ➔ Make sure that the bed stands in straight position before each weighing operation. Check this if necessary with the circular level.
- ➔ Initial calibration of the weighing system is performed in the factory. No recalibration is required for the intended application environments 3 and 5 (see chapter 2.3).

6.6 Weight monitoring

- 01** The saved weight change values can be displayed as a curve on the CareBoard® display.
- 02** Save & Send: Save and send the weight currently displayed (transmission of weight data to the in-house care software).
- 03** The left field on the display shows the absolute weight in steps of 0.5 kg.
- 04** The right field on the display shows the increase and loss of weight in steps of 0.5 kg.
- 05** Press this button to enter the resident's height. Based on the weight measured, the resident's BMI is then calculated automatically and displayed next to the button.
- 06** With the freeze button, you can add additional weights, or remove weight from the bed without influencing the absolute weight or the weight change.
- 07** Here, you can set the absolute weight and the weight changes to 0 kg.
- Press once and the right field on the display flashes.
 - Press a second time and the value is set to 0 kg.
- If the resident is in bed:
- Short pressing: The relative weight display flashes
 - Pressing once again: Relative weight display is set to zero
- In case of a change of resident:
- Long pressing: Both displays flash
 - Pressing once again: Both displays are set to zero


08 By using the +/- button, you can make the right display of the CareBoard® show the difference between the weight currently displayed and the saved weight. After selecting this panel, the display is visible for approximately 5 seconds.

- ➔ If installed at the foot end, the Careboard® must be taken out of its holder for weight monitoring to prevent falsified display values.
- ➔ To ensure accurate weighing, the weighing system is to zero before use. For this purpose, the bed must be unoccupied.

If you have to add or remove additional weight such as cushions or other items but nevertheless want the resident's actual absolute weight and weight change to stay unaltered, follow the instructions described below:

- Press the freeze button.
Both displays are showing the freeze symbol.
- Now put the additional weight on the bed or remove it, and press the freeze button once again.
The absolute weight and the weight change are shown unchanged again.

6.7 Weight history

01 Return to weighing menu.

02 Here you can select the displayed time period: one day, one week or one month.

03 Here you can delete all saved values of the weighing system in case of a resident change.


04 Here you can open the saved values of the previous period of time.

- In order to delete the weighing system's values, press the delete button for approximately 5 seconds.

During this time, a white bar next to the button gets smaller.

6.8 Nursing routine documentation

The nursing routine documentation allows resident data to be entered directly by the bed. It also comprises a checklist for the daily nursing routine. The Save & Send function automatically transmits the data to the in-house care software. This is only possible if the nursing software used supports the data import and the home's own WLAN/LAN network allows transmission.


01 Forward/Back

02 Save/send data

03 Keypad for data input

04 Vital data display

- ➔ Pressing the floppy disk starts the saving/transmission process. Meanwhile it flashes yellow. The data then disappears from the fields and can be viewed in the in-house care software. If the send process fails, an error message will appear and the Save & Send button flashes red.


01 Check box for daily nursing routine

01

6.9 Nurse settings

This menu can be used to set the time, the language, and the call volume (at the CareBoard®) as well as the volume for the keypad tones.

01 Activate or deactivate the clock (screen saver)

02 Set the date and time

03 Select language

04 Delete resident data


05 Set the volume of the nurse call (of the CareBoard®)

06 Display the current date and time

07 Set the button sound volume

6.10 Underbed light (UBL)

The underbed light can be permanently switched on, off or can be controlled automatically. When automatic control is activated, the underbed light switches on whenever the resident sits on or near the edge of the bed and leaves the bed. It switches off when the bed is occupied again. Please activate the automatic function for when the resident is lying in bed. In that moment the function is calibrated.

01 Switching on the underbed light

02 Switching off the underbed light

03 Automatically switch on/off the underbed light


If the resident rests near the edge of the patient surface, the automatic function could see this as a preparation to get up and switch on the light. If this happens frequently and the resident is bothered by this, the automatic function should not be used.

6.11 Digital user manual

In addition to the printed instruction manual, you can also access a digital user manual when using the CareBoard®. It provides you with an explanation of the current screen.

6.12 Service reminder (pop up)

The CareBoard® reminds the user of the service date two weeks before it is due. This message is displayed until it is actively acknowledged by pressing „OK“. On the service date, the message „Service overdue“ is displayed, which also has to be acknowledged by pressing „OK“.

➔ If you do not wish to get a service reminder, you can have it deactivated in the CareBoard® by a service engineer.

7. Accessories

7.1 Mattress

The bed is equipped with a mattress. If necessary, this mattress can be replaced by the same or a similar product.

Caution!

Do not use water mattresses. If the safe working load is exceeded, the drives and other components may fail and cause accidents. Never exceed the safe working load of 250 kg.


Caution!

Only replace wissner-bosserhoff mattresses with equivalent mattresses of the same dimensions. Make sure that only mattresses with incisions on the underside that allow for movement are used, as a mattress without these could affect the adjustment functions. Do not use spring mattresses.

We recommend to use mattresses manufactured with a minimum volume weight of 40 kg/m³ (RG 40). All wissner-bosserhoff standard mattresses fulfill this requirement. Please refer to the following overview for the suitable mattress dimensions for the corresponding patient surfaces.

Patient surface type	Dimensions
Aero / Easyclean	12-18 cm × 87 cm × 200 cm
Comfort	12-15 cm × 87 cm × 200 cm


7.2 Lifting pole

The holder for the lifting pole is located at the head end of the patient surface frame.

Attention!

Grab handle and belt consist of high-quality plastic. However, all plastics are not unconditionally resistant to aging, so that grab handle and belt can break or tear after long usage. This can cause injuries. Replace grab handle and belt after 4 - 5 years at the latest. The year of manufacture is imprinted on the product.

If you want to attach a lifting pole:

- From above, insert the lifting pole into the holder. It is held in its high position by a bolt.
- Turn the lifting pole, so that the bolt lies in one of the notches provided in the holder. Vertical and horizontal adjustment is possible.

Attention!


Under high load, the grab handle and belt can break or tear and cause injuries. Only use a patient lifting pole under a maximum load of 750 N (approx. 75 kg).

- Attach the grab handle to the patient lifting pole, using the belt.

7.3 Infusion stand

The infusion stand is only intended for attaching infusions. Do not attach infusion pumps here.

→ The infusion stand is only intended for attaching infusions. Do not attach infusion pumps here.

The holder for the infusion stand is located at the head end of the patient surface frame.


If you want to attach an infusion stand:

- From above, insert the infusion stand into the holder.
- Turn the infusion stand, so that the bolt lies in one of the notches provided in the holder.

7.4 Extensions of patient surface

The low nursing bed sentida 7-i offers two methods to extend the patient surface. A permanent extension of the patient surface, carried out at the manufacturer's plant, or an optional temporary extension.

Long-term (fixed) extension for the patient surface


Can be folded down

The low nursing bed sentida 7-i can be delivered with a 10 cm or 20 cm extension for the patient surface. This patient surface extension can also be implemented afterwards in the factory or by authorized qualified personnel. An additional extension can then be attached to the calf rest. To close the resulting foot-end gap, a corresponding protector needs to be inserted. If you want to have the patient surface of your low nursing bed sentida 7-i extended, please contact your service partner. You can find the address at the end of these instructions for use.

Temporary extension of the patient surface

If the bed is equipped accordingly, the patient surface of the bed can be extended and shortened again by 10 cm, 20 cm or 30 cm at the foot end without tools. This is why the foot end of the bed is retractable and can be locked at two levels. The foot end can be unlocked with a yellow rope, which is attached between the side panels at the left and right lock of the foot end. You can unlock the foot end by pulling the yellow rope slightly upwards. You need to insert a protector to close the resulting foot-end gap in this case as well.

Foot-end

Yellow rope for releasing
the foot end

**Danger!**

If you fail to close the foot-end gap resulting from the patient surface extension with a protector, there is danger of entrapment! This applies to both the temporary and the permanent patient surface extension.

If you want to extend the patient surface of the low nursing bed sentida 7-i:


- Pull the yellow rope at the inner side of the foot end slightly upwards.
The foot end is now unlocked.
- Pull on the foot end until you hear a clicking sound at the desired position with the rope being loose.

In order to shorten the low nursing bed sentida 7-i again, just follow the steps in reverse order:

- Unlock the foot end by pulling the yellow rope and push it back in.ein.

7.5 Foldable linen holder

At the outer side of the foot end, a linen holder can be attached as accessory. The linen holder is connected to the foot end with a folding hinge. It serves for temporarily putting linen or other objects down during the nursing care process. Assembly instructions are available from wissner-bosserhoff for later extension.

- To open the linen holder, grab it at the top, pull it upwards and flap it down.
- To close the linen holder, fold it upwards and lower it in the folding hinge.

**Attention!**

With moving accessories, there is a risk of entrapment. Move the accessories carefully in order to avoid squeezing.

Do not load the linen holder with more than 15 kg.


Please note that the linen holder cannot be used simultaneously to the CareBoard® being attached at the food end of the bed.

7.6 Integrated linen holder

The integrated linen holder is permanently installed in the lying surface and can be pulled out through a recess in the footboard. It serves for temporarily putting linen or other objects down during the nursing care process. Assembly instructions are available from wissner-bosserhoff for later extension.

- To pull out the linen holder, hold the aluminum bar and pull as required. If you need it, you may fold out the bracket.
- To retract the linen holder, fold down the bracket and carefully push the linen holder in until it clicks into place.

Attention!

With moving accessories, there is a risk of entrapment. Move the accessories carefully in order to avoid squeezing.

Do not load the linen holder with more than 15 kg.

7.7 Accessories holder

The accessories holder serves for hanging urine bag holders, urine bottle holders or similar accessories on it. It is attached to the side panel.

The maximum load on the accessory holder is 2 kg.

7.8 Telescopic wall spacers

The telescopic wall spacers are available in short, medium, and long versions for both under-carriage types (homely and extra mobile).

head side

lateral

The wall spacers are attached to the undercarriage of the bed with screws.

- Open and remove the tube clip to adjust the length, pull out the telescopic arm to the desired length, and then reinsert and tighten the tube clip.

7.9 Side guard padding, side guard net

Protective padding or netting are available for the SafeFree® side guard. These can provide the resident with additional protection by reducing or closing of gaps.

- To attach the padding, place it against the side guard element and close the snap buttons at the bottom of the padding.
- To attach the netting, place it against the side guard element and feed and attach the Velcro fastener tabs around the side panel.


7.10 Extension for SafeFree® side guard (accessory)

In order to increase the height by approx. 60 mm, an extension is available as accessory for the SafeFree® side guard.

Attention!

An unsuitable extension, which is not intended for the SafeFree® side guards, could come loose and cause a fall of the patient or resident. Also ascertain yourself that the extension is suitable for the SafeFree® side guard. Ensure the secure attachment of the extension on the SafeFree® side guard.


If you want to attach an extension to the SafeFree® side guard,

- Completely lift the SafeFree® side guard.
- Attach the extension at the center of the side guard's handrail.
- Press down the mounting clamps on the handrail until they engage and fix the extension.

If you want to remove an extension from the SafeFree® side guard,

- Position both SafeFree® side guards at different heights.
- Press the SafeFree® side guard's release button and push the extension out over the hand rail.

8. Cleaning and disinfection

The nursing home bed must be cleaned and disinfected at regular intervals (e.g. as specified in the hygiene plan), as circumstances dictate (e.g. soiling after a particular episode) or at least before each change of resident as per the hygiene plan.

Danger!

If electrical components come into contact with fluids, this can significantly increase the risk of electric shock or fire. This means you must disconnect the bed from the power supply prior to cleaning/disinfection.

Attention!

Damaged surfaces should be repaired/replaced immediately as they cannot be cleaned/disinfected properly and allow fluids to penetrate. This creates a risk of infection and damage to the bed.

Please note!

The nursing home bed may not be cleaned/disinfected in automatic washing units or with spray water. Only wipe cleaning/wipe disinfection is permitted.

Please note!

If you are unsure, check first whether the cleaning agent/disinfectant is suitable for the relevant surface by applying some to an area out of sight.

Please note!

Following cleaning/disinfection, it is important to ensure things are properly neutralized and there is no moisture or cleaning agent residue on any surfaces.

8.1 General procedure for cleaning/disinfection

- Disconnect the bed from the power supply.
- Make sure the product is free from any obvious dirt or particles of dirt that might damage the surface during cleaning/disinfection.
- Wipe the bed thoroughly with a soft cloth that has been lightly moistened with the all-purpose cleaner/disinfectant solution.
- Then wipe the product thoroughly, to neutralize it, with a soft cloth that has been lightly moistened (with water only)
- Make sure the product is free of any cleaning agent or disinfectant residue.
- Dry the product with an absorbent, lint-free cloth.
- Make sure there is no moisture/humidity left on the surfaces or in any gaps.
- Reconnect the bed to the power supply.

8.2 Cleaning

The surfaces of wissner-bosserhoff furniture products comply with the requirements of EN 12720 for the surfaces of furniture in terms of resistance to chemicals and are therefore highly resilient to the normal kinds of stress they might typically encounter.

It is essential, however, not to use any cleaning products or cleaning agents that contain corrosive or abrasive ingredients and act like a scouring agent.

Only all-purpose, weakly alkaline cleaning products (that create soap suds), with tensides or phosphates as the active cleaning ingredient, should be used. Dose the all-purpose cleaner in accordance with the specifications of the cleaning agent manufacturer.

8.3 Disinfection

Disinfection will only be effective if the bed has been cleaned beforehand.

Only suitable disinfectants equivalent to the disinfectants described in Annex 1 to EN 12720 Part 1 or the agents in the list drawn up by the German Society of Hygiene and Microbiology or the Robert Koch Institute may be used. The disinfectant should be used (e.g. concentration) precisely as specified by the respective disinfectant manufacturer.

If you have any doubts about the suitability of a given disinfectant, please contact us.

8.4 Special information regarding real wood parts

We have used materials of the highest quality to create a product for you with a high level of benefits and homely real wood character. Slight differences in shade, contrast and color between veneer surfaces, solid wood parts or decorative surfaces (melamine, HPL), as well as anything showing up in real wood materials as a result of the growth process, do not count as defects.

Varnished real wood parts (profiles at the head and foot ends) need to be cleaned/disinfected with special care as they are very sensitive to moisture. No moisture may be left on these surfaces or this may damage the varnish and wood.

9. Maintenance and service

In accordance with

- the applicable national regulations for medical devices and medical beds
- the current national health and safety regulations and national repetitive testing standards for electrically operated equipment

operators of beds for medical use are required to ensure the safe condition of the medical devices throughout the entire service life. This includes professional maintenance and safety inspections at regular intervals.

The nursing home beds from wissner-bosserhoff GmbH are designed for a service life of 10 to 15 years. The beds require very little servicing. This is because, during product development, care was taken to reduce that servicing work to a minimum, with low operating costs.

In daily use, however, experience shows that products are also sometimes handled carelessly, and rough handling can also contribute to them aging more quickly and to the wear of certain components, on which the manufacturer can have no direct influence.

The manufacturer assumes liability for the safety and reliability of the product only if it is regularly maintained and used in accordance with the operating, usage and safety notes in these instructions for use.

If severe defects that cannot be rectified are detected in visual inspection, functional testing, checks, measurement or during maintenance, the product must be taken out of use. In addition to the comprehensive inspections by professional technicians, users (nursing staff, family carers, etc.) should perform visual inspections and functional tests at shorter, regular intervals. This is especially important when a bed is used for a new resident.

RECOMMENDATION: The safety inspection in accordance with IEC 62353 should take place every 12 months in conjunction with mechanical maintenance.

Observe the following order of inspections and tests according to IEC 62353-VDE 0751-1:

1. Visual inspection
2. Electrical measurement
3. Functional test

Visual inspection and functional test

- In accordance with the IEC 62353-VDE 0751-1, these inspections/tests must only be performed by qualified persons who have received adequate training and have access to the required equipment.

Electrical measurement

- Electrical measurements according to IEC 62353-VDE 0751-1 may also be performed by specially trained and technically authorized persons (e.g. in-house technicians) with additional medical and bed-specific knowledge if suitable measurement equipment is available.
- Only a qualified electrician with additional medical and bed-specific knowledge may evaluate and document the test results.

IEC 62353-VDE 0751-1 a short checklist

Visual inspection

- Mechanical components: (undercarriage, side guards, wooden frame, patient surface frame, lifting pole socket)
- Electrical components: power and handset cable, handset
- Leakage current test (difference ΔI without load)
- Protective conductor testing (protection class 1)

Functional test:

- Mechanical components: castors, emergency lowering of the backrest, side guards, calf rest, motor bolts, screws
- Electrical components: Handset

Documentation:

- Test protocol, non-conformity report, inventory list, instructions for use

9.1 Spare parts

Danger!

There is a highly increased risk of electric shock or fire when using incompatible electric components. Only replace electrical components with original spare parts.

Only wissner-bosserhoff GmbH original spare parts may be used. Contact our customer service, sales and technical consultation department for information concerning spare part deliveries. On request, please provide the details on the nameplate and the order number.

Our qualified technical customer service will be happy to help:

Phone +49 2377 784-456
Fax +49 2377 784-150

9.2 Troubleshooting

Malfunction	Possible cause	Remedy
Electric motor adjustments do not function correctly.	Mechanics are locked.	Check moving parts and remove any foreign bodies.
Electric motor adjustments do not function.	Handset faulty or not plugged in correctly.	Contact our customer service and request a replacement.
	Plug not connected to the power supply and battery (if applicable) empty	Connect the plug with the power supply system.
	The system was overloaded and the temperature protection has switched off; the fuse has tripped	Wait for the system to cool down (up to 20 min.). If the system still does not function, contact our customer service. (the power box may need to be replaced)
The SafeFree side guard cannot be adjusted properly anymore.	Mechanics are locked.	Check moving parts and remove any foreign bodies.
	Mechanical parts are bent	Please contact our customer service.
The castors do not brake or do not roll.	Impurities have become caught in castors over time.	Remove impurities.
	The castor system is defective.	Please contact our customer service.
Bed does not move/moves in an undesired way.	Software or hardware of the CareBoard® is defect	Reset the bed by pushing the STOP button; in case of a permanent error, please contact the WiBo customer service
Alarm is not triggered	Defect sensors for humidity, brakes and side guards	Contact the customer service
	Humidity sensor mat not connected correctly	Check the correct connection of the mat
	The sensitivity of the sensors is set incorrectly	Make new settings while the resident is lying in bed
Weighing system menu is not displayed correctly	Weighing box is defective or not connected correctly	If it is suspected that the weighing box is defect, please contact the WiBo customer service.

Alarm is triggered without any apparent reason	Nursecall system adapter is defective	Contact in-house technician, otherwise contact the WiBo customer service
	Cable is defective	
	Sensors are defective	
	Missing magnets	
Weight is not displayed correctly or values of the weight history differ too much	Weight was added/removed without freezing it	Ignore extremely differing/ wrong values
	Weight was set to zero accidentally	
	Wrong measuring	
Display shows „hi“ or „lo“	Zero point set incorrectly	Set the display to zero when the resident is not in bed
Underbed light is not working as it should	Wrong standard saved (patient is in the bed when turning on automatic functions)	Switching underbed light into automatic mode when the resident is in bed
Handset does not work	Handset locked, not correctly plugged in or faulty	Unlock handset, check connection. Otherwise contact the customer service
Mobilift® does not work	GO button is not pressed	Press the GO button (handset or CareBoard®)
	Handset is locked	Unlock handset
	Mobilift is defective (cable break)	Contact the WiBo customer service
CareBoard® does not react properly	You are pressing the touchscreen too hard or the area being pressed is too large	When operating the board, please make sure that you press the center of the button and do not push too hard
The power box emits an acoustic signal	CareBoard® not correctly connected / faulty	Check the plug connection or call WiBo service
Save & Send button lights up red	Fault affecting the server connection	Check network cable, reset the CMC box (see chapter 5.2) – or contact WiBo service

9.3 Storage

When stored, the bed may be exposed to environmental conditions that are within the following limit values:

Ambient temperature	+10 °C to +40 °C
Relative humidity	30 % to 75 %
Air pressure	800 hPa to 1060 hPa

Caution!

If the sentida low nursing bed has been stored or transported at low temperatures, it needs a certain time and temperature to adjust before use. If the adjustment time is too short or the temperature is not suitable, the sentida low nursing bed can be damaged and fail. After strong temperature fluctuations, let the sentida low nursing bed adjust its temperature for at least 12 hours.

For storage of the low nursing, make the following preparations:

Unplug the power supply.

- Remove all accessories such as the lifting pole, serving tray etc.
- Remove the battery.
- Cover the low nursing bed so that it is protected from damage.
- Note down the storage date so that you can comply with the maintenance intervals.

10. Disposal

This device falls within the scope of the EC Directive 2002/96/EC (WEEE). It is not registered for use in private households. Disposal at municipal collection points for discarded electrical devices is not permitted. wissner-bosserhoff GmbH is responsible for legal compliance concerning the disposal of this device. For further information, please contact your responsible sales partner, or our company, if you are situated within Germany.

When transferring the equipment to an industrial third party, you are bound by contract to point out that after the end of use, proper disposal must be undertaken or arranged. If you do not point this out to the third party, you are responsible for the proper disposal of the equipment after the end of use.

The low nursing bed contains lead gel batteries, electrical parts and metal parts. It can also contain plastic parts made of ABS, PA, PUR or PE.

The metal and plastic parts that accumulate during service and repairs must be disposed of properly and professionally in accordance with the applicable laws and regulations. Specifically the electric motors and electrical controls may only be disposed of through specialist firms authorized to do this, or through the respective waste management facilities.

11. Warranty

Statutory warranty provisions apply.

This warranty covers all malfunctions and defects associated with materials and manufacture. Malfunctions and defects arising from incorrect handling and external influences are excluded from this. If there are grounds for justified complaints during the warranty period, these shall be dealt with free of charge. This warranty can be enforced with the sales receipt, which bears the date of purchase. Our general terms and conditions of delivery apply.

12. Technical data

Outer dimensions	approx. 208 x 106 cm (depending on the wooden lining)
Nominal dimensions of patient surface	90 cm × 200 cm
Height adjustment	27.5 cm – 80 cm (measured at the seat section)
Backrest adjustment	65° ± 5° incl. backrest retraction by approx. 11.5 cm
Thigh rest adjustment	30° ± 5° incl. thigh rest retraction by approx. 4.8 cm
Tilting	Head-end approx.15° / foot-end approx.15°
Safe working load	250 kg (215 kg resident with 35 kg mattress + accessories) *
Bed weight	ca. 145 kg (depending on configuration)
Type of protection	IPX4
Protection class	II
Operating time	INT 2 min./18 min. The equipment is not intended for continuous operation. After 2 minutes of operation at full load,a break of 18 minutes is recommended.
Power cable	CAT 5
Noise level, maximum	ca. 57 dB/A
Application environment	3 und 5
Relative humidity	30% – 75%
Atmospheric pressure	800 hPa – 1060 hPa
Environment temperature	+10°C – +40°C
Product life	If used, cleaned, serviced and repaired correctly, a useful life for our products of 10-15 years can be expected. Excluded from this are wear parts such as castors, electrical components, etc.
Manufacturer	wissner-bosserhoff GmbH Hauptstraße 4 – 6 58739 Wickede (Ruhr) GERMANY Tel. +49 2377 784-0

* **Manufacturer's recommendation:** the maximum resident weight may be able to be increased, depending on the weight of the mattress and accessories used. For example: safe working load 250 kg - mattress 10 kg - accessories 10 kg = max. resident weight 230 kg. When using heavier mattresses/ accessories, the max. permissible resident weight changes accordingly!

Power-Box	
Electrical connection	230 V/50 Hz
Type of protection	IPX4
Protection class	II
Operating time	INT 2 min./18 min.: The equipment is not intended for continuous operation. After 2 minutes of operation at full load, a break of 18 minutes is recommended.
Power consumption	up to approx. 280W

Height adjustment motors	
Electrical connection	24 V DC
Type of protection	IPX4
Protection class	III
Operating time	INT 2 min./18 min.: The equipment is not intended for continuous operation. After 2 minutes of operation at full load, a break of 18 minutes is recommended.
Adjustment speed	approx. 5.4 mm/sec.
Pressure force	max. 6000 N

Patient surface motors	
Electrical connection	24 V DC
Type of protection	IPX4
Protection class	III
Operating time	INT 2 min./18 min.: The equipment is not intended for continuous operation. After 2 minutes of operation at full load, a break of 18 minutes is recommended.
Adjustment speed	10 – 14 mm/s
Pressure force	max. 3000 N - backrest motor with emergency lowering max. 6000 N - backrest motor without emergency lowering (optional) max. 6000 N - thigh rest motor

CMC-Box	
WLAN connectivity	802.11 b/g/n
Frequency	2,4 GHz
Encryption	WPA2
Ethernet port	10/100 MBit/s
Protocols	IPr4 TCP

Radio module RTM08	
Frequency	868,3 MHz
Transmission power	0,25 mW
Supply voltage	5 V DC
Standby current	max. 1,5 mA
Transmission current	max. 13 mA

13. Classification

The low nursing bed sentida 7-i is a medical product class Im (non-sterile with measuring function) in accordance with appendix IX of the directive for medical devices 93/42 EEC and amending directive 2007/47/EC.

13.1 Applied standards

Designation	Comment
Directive 93/42/EEC	Medical devices directive 93/42/EEC
German law on medical products	German law on medical products (national implementation)
DIN EN ISO 14971:2013-04	Risk management applied to medical products
DIN EN 60601-1:2006 + A1:2013 (applicable sections)	Medical electrical appliances
DIN EN 60601-2-52:2010 + A1:2015 (relevant sections)	Medical beds
BfArM [German Federal Institute for Drugs and Medical Devices] recommendations	Recommendations of German Federal Institute for Drugs and Medical Devices
DIN 33402-2:2005-12 + amendment 1	Human body dimensions
DIN 68861-12011-01 DIN EN 12720:2009-07	Furniture surfaces subjected to chemicals, resistance to cold liquids.
DIN EN 60529; VDE 0470-1:2000-09	Type of protection by housing IP code (protection from humidity)
DIN EN 12530:1999-05 (relevant sections)	Castors - equipment castors - manually controlled castors
EN 50419:2006/ directive 2002/96/EG (WEEE)	Labeling of electric and electronic equipment
Directive 2014/31/EG	Directive on non-automatic scales
EN 45501 : 2015	Metrological aspects of non-automatic scales
EN ISO 15223-1	Symbols to be used with medical device labels


14. Product identification

14.1 Symbols used


	General warning!	Black exclamation mark in yellow triangle
	Follow the instructions for use!	White person on blue background
	Caution: Possible cable entrapment hazard!	Black power cable in yellow triangle
	Caution: Possible foot entrapment hazard!	Black foot in yellow triangle
	Caution: Possible hand entrapment hazard!	Black hand in yellow triangle
	Type B applied part	
	Das Gerät unterliegt unter anderem der EG-Richtlinie 2002/96/EG (WEEE). Es wurde nach dem 13.08.2005 in den Verkehr gebracht.	
	Conformity mark for the directive on scales	Black „M“ on green background
	Earth (functional earth)	
	Precision class of scale	
	Radio symbol	
	Adult detection, minimum weight, minimum height, minimum BMI	

14.2 Nameplate

1	Manufacturer's logo
2	Manufacturer's address and Homepage
3	Model name
4	Article number
5	Serial number
6	Safe working load
7	Operating time (2 min. operation, 18 min. pause)
8	Input voltage/frequency/apparent power
9	Mandatory sign „Observe instructions for use“
10	Codes for protection rating in accordance with DIN 40050: 1. code: Protection rating for contact and foreign body protection 2. code: Protection rating for protection from water (4= spray water)
11	Radio Symbol
12	Signal grounding (for electromagnetic interference suppression, no protective conductor function)
13	Protection class II device (classified according to the type of protection from electric shock)
14	Caution: Electrical waste - statutory requirements for the disposal of old equipment must be complied with!
15	Type B applied part (classified according to the degree of protection from electric shock)
16	For indoor use only
17	CE mark including year of marking and identification numbers of notified bodies for conformity assessment

14.3 Scale nameplate

1	Factory brand, name and address of manufaturer
2	Precision class of scale (III = commercial scale)
3	Number of calibrations
4	Serial number
5	EC design approval
6	Calibration value (in compliance with DIN EN 45501)
7	Tare maximum load
8	Minimum load
9	Maximum load
10	Permissible operating temperature
11	Warranty seal
12	Type designation of scale

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Reproduction, also excerpts, only with prior permission of the publisher. All rights reserved.