SafeSense®

Bed-Exit System

Gebrauchsanweisung
Instruction manual
Notice d'utilisation
Gebruiksaanwijzing
Istruzioni d'uso
Manual de Instrucciones

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1. Table of contents

1. Table of contents	30
2. Introduction	32
2.1 About these instructions for use	32
2.2 Symbols used in the text.....	32
2.3 Intended use.....	33
2.4 General safety notes	34
3. Components.....	35
4. System structure.....	36
4.1 Cable version	36
4.2 Wireless version	36
5. General usage information	37
6. Start-up.....	38
6.1 Delivery	38
6.2 Installing on a nursing home bed with motor „Duomat 9“	38
6.3 Installation on a nursing home bed with motor „Quadromat 7“	40
6.4 Installing on a nursing home bed with motor „Quadromat 7“ and battery	41
6.5 Cable routing for nursing home beds	42
6.6 Securing the sensor mat.....	43
6.7 Connecting to the nurse call system or other signal receivers	44
7. Operation of bed exit box.....	45
7.1 Switching on/off.....	45
7.2 Setting a time interval	45
7.3 Display	45
7.4 Underbed light	46
8. Permissible mattresses and minimum resident weight.....	47
9. Cleaning the sensor mat	48
10. Maintenance and service	48
11. Spare parts	49
12. Disposal	50
13. Warranty.....	50
14. Technical data.....	51

15. Classification.....	51
17. Nameplate	56
18. Transport and storage.....	56
19. Troubleshooting.....	57
20. Declarations of Conformity.....	170

2. Introduction

2.1 About these instructions for use

In this section, you will find information about the structure of these instructions for use and a description of the signs and symbols used.

This document contains instructions for the use of the SafeSense® Bed-Exit System.

These instructions for use may contain inaccuracies and typographical errors. The information provided in these instructions for use are updated at regular intervals, and changes resulting from product maintenance are implemented in future editions. Changes or improvements are possible at any time without previous notice. Contact our customer service in case you have any questions.

The instructions for use must be read and observed by every person operating the Bed-Exit System.

Apart from the instructions for use and the accident prevention regulations valid for the respective country and area of use, the commonly accepted regulations for safe and professional work must also be adhered to.

2.2 Symbols used in the text

In these instructions for use, the following terms and symbols are used for important notes:

Danger!

Safety notes on people being in danger are marked with this symbol. This symbol is used when there is imminent danger of death or severe injuries.

Attention!

This symbol is used when there might be the danger of unsafe situations or minor injuries.

Caution!

This symbol precedes warnings when there is the danger of damage to the equipment or other things.

This symbol precedes additional helpful advice.

- A hyphen preceding the text means: This is part of an itemization.
- A dot preceding the text means: You have to adhere to this instruction.
Indented text means that this is the result of your actions.

2.3 Intended use

SafeSense® is classified as a Class I medical device as an accessory for medically used beds. It is intended exclusively for human medical and nursing purposes and for use in the physical environment of nursing homes, outpatient care facilities and clinics in accordance with application environment 3 and/or 5 as per IEC 60601-2-52 (current status). SafeSense® may only be used in combination with the wissner-bosserhoff nursing beds of the sentida and movita model series and hospital beds of the image 3 and eleganza 2 model series.

We assume no liability for potential property damage or personal injuries caused by unapproved accessories, inappropriate handling, accidental tripping or mutual cancellation of the intended use. The specified use is the intended use. Its sole source for the operator or user is the labeling and the instructions for use.

Indication and target patient group

SafeSense® is intended for monitoring general presence parameters in bed. It informs the nursing staff about leaving the bed. At night, it provides the resident with automatic underbed-light.

SafeSense® is intended for persons who are at least 146 cm tall, weigh at least 40 kg, and have a BMI ≥ 17 . Reliable signal transmission can only be guaranteed under these conditions.

Other applications must be agreed in writing with wissner-bosserhoff GmbH in advance. The product is to be used as work equipment for care applications and is subject to the regulations of the responsible professional associations.

Application restrictions and contraindications

Use only under medical supervision. The determining factor as to whether or not care is medically supervised is whether that care is provided under the direction of medical personnel. SafeSense® is not intended for use as an emergency call or lifesaving call for help. Rather, it is a tool to facilitate everyday care activities in nursing homes.

Intended users

The system may only be used by persons who, based on their training or knowledge and practical experience, can guarantee proper handling. In addition, the user of the system must have been instructed in its correct use and must have familiarized himself with the product on the basis of these instructions for use. Common users include: nursing staff, physicians, physical and occupational therapists and house technicians.

2.4 General safety notes

The SafeSense® bed exit system was designed in accordance with state-of-the-art technology and approved safety regulations.

Use SafeSense® only when it is in perfect working order. Only use it in for its intended purpose, safely, risk-consciously and in accordance with these instructions for use. Malfunctions that may affect safety must be eliminated immediately in particular!

Always keep these instructions for use ready to hand at the location where the Bed-Exit System is used. In addition to the instructions for use, please observe the generally applicable standards and commonly accepted regulations on accident prevention and environmental protection!

Do not carry out any changes, extensions or conversions without the manufacturer's approval. Spare parts have to comply with the manufacturer's requirements. For original spare parts, this is always guaranteed.

Make sure that operating and auxiliary materials and replacement parts are recycled safely and in an environment-friendly manner.

3. Components

The SafeSense® Bed-Exit System is available in a cable version or a wireless version and consists of the following components:

Bed-Exit Box

Sensor mat, 80 x 20 cm


Wireless receiver for nurse call system
(wireless version only)


Customer-specific adapter cable with
connection to the nurse call system (with
optional on/off switch)

4. System structure

4.1 Cable version

4.2 Wireless version

5. General usage information

When using SafeSense®, please be aware that there are various circumstances that could lead to malfunctions or complete loss of function in the system. Therefore, please note the following before and during use:

- Check that all plugs are connected correctly to ensure that all signals can be transmitted.
- Make sure the system is switched on
- Ensure that all cables are undamaged
- Position and fasten the sensor mat correctly to keep it from slipping (see “Securing the sensor mat”)
- It must always be ensured that the resident/patient is positioned correctly in the bed (e.g. not facing backwards)
- Please observe the relevant minimum weights associated with the various mattress types (see “Permissible mattresses and minimum resident weight”)
- If incontinent residents/patients occupy the bed, incontinence protection must be used on the mattress
- Do not use the sensor mat in combination with alternating pressure mattresses. Otherwise, the system will not be able to transmit reliable signals
- Consider that the sensor mats are also subject to wear and tear to the degree that can be reasonably expected from their use and must be replaced when they begin to show signs of functional defects.*
- Keep the sender and receiver together if possible. If this is not possible, delete the receiver's programming (the exact steps can be found in its usage instructions) and reconfigure it

To identify possible errors, please also note the error-cause matrix given later in the instructions for use (see chapter “Troubleshooting”). This lists faults that may occur along with their possible causes.

* The wear and tear depends on the resident's/patient's weight, the bed occupation time and the strength and other properties of the mattress. The sensor mat has to be replaced after 5 years at the latest.

6. Start-up

6.1 Delivery

The SafeSense® Bed-Exit System is generally delivered as individual components and is assembled by specialists on site.


- Please check whether all components have been delivered by consulting the delivery papers.
- Write down potential defects or damage on the delivery note.
- Inform your service partner immediately about potential transport damage or defects. The address and phone number are noted on the last page of this manual.

6.2 Installing on a nursing home bed with motor „Duomat 9“

Use the handset to position the bed at a comfortable working height and pull the power plug.

Attach the retaining plate to the bed exit box and tighten the screws. The wide holding tab should point down towards the underbed light.

Use a screwdriver to remove the protective motor cap. (There is one screw on the left side and one on the right)


Disconnect the handset (third plug from the left) and connect the bed exit box.

Open the fastening clip of the motor and slightly pull out the slider. Attach the box to the motor by inserting the holding tab from below, hook the upper holding tabs into the slider and put it in the initial position. Close the fastening clip.

Then, put the protective cap back onto the motor.

➔ **Fasten the box to the motor under the backrest at the head end of the bed!**

Plug the handset and sensor mat into the bed exit box (the handset on the left, the sensor mat on the right). Attach the strain relief for the handset cable.

Insert the nurse call cable in the bed's cable conduit and fasten it to the strain relief fixture with a cable tie. Reconnect the power plug.

➔ **Scan this QR code to see the step-by-step installation video:**

6.3 Installation on a nursing home bed with motor „Quadromat 7“

The SafeSense® system is quick and easy to install. To do so, please follow the steps below:

Disconnect the bed from power supply. Use a screwdriver to remove the drive cap. There is one screw on the left side and one on the right side. Remove the handset from the third socket from the right.


Attach the retaining plate for the motor and tighten the screws (position the plate so that the holding tabs are on the top and on the bottom)

Connect the hand switch to the bed-exit box and screw on the ignition relief as shown in the illustration.

Clip the Bed-Exit Box into the holder. The display has to point upwards when you do this. The box will sit securely once it has audibly slotted into place. For a bed with a battery, the box is screwed to the battery and slotted together with it. (see chapter 6.3)

Connect the Bed-Exit box to the bed drive. To do this, plug the connector plug (on the left of the box) into the now empty slot on the drive.

Insert the sensor mat into the Bed-Exit Box. To do this, first remove the plug of the connection socket below the dial for call interval setting. Then position the sensor mat as described in chapter 6.6. Reconnect the power plug.

For the wireless version (see photo), please plug the wireless receiver's power supply into a power socket and connect the plug to your call system. The wireless receiver must then be calibrated to the SafeSense® system.

The final step is to connect the SafeSense® system to your nurse call system. For the cable version, this is done by connecting the plug on the SafeSense® system to the socket on your call system.

6.4 Installing on a nursing home bed with motor „Quadromat 7“ and battery

1. Remove the battery pack from the bed and remove the two diagonally opposite screws.

2. Attach the retaining plate to the bed exit box. (The holding tabs should be on the sides)

Secure the Bed-Exit Box to the battery using the two screws supplied.

Clip the battery into the holder on the bed.

The remaining installation steps are identical to those for beds with no battery.

6.5 Cable routing for nursing home beds

When installing SafeSense® please take care to ensure that the cables are not pinched or damaged in any way. For this reason, the steps below are to be followed to ensure proper cable routing:

1. Cable routing for beds with cable opening on the side:

Use a cable clip to attach the nurse call cable to the power cable.

Always run the nurse call cable along the power cable to the head end of the bed and attach it to the power cable using a cable clip.

2. Cable routing for beds with cable opening at the head end:


Run the nurse call cable from the Bed-Exit Box into the cable conduit, which also contains the power cable. Use a cable clip to attach it to both the inside of the head end and the cable relief at the head end.

3. Cable routing for beds with a D-Box on the frame:

Guide the nurse call cable to the D-box and insert it into the cable conduit.

4. Cable routing for beds with a D-Box at the head end:

Run the nurse call cable along the power cable and attach it to the holder in front of the cable conduit, which contains the power cable. Then run the nurse call cable through the cable conduit and below or above the D-Box. Use a cable clip to attach the cable. Continue running it alongside the power cable and tie it to the power cable at the point where it starts to coil.


Attention!

Please note: Gather the remaining cable together to avoid the risk of tripping and cable damage.


6.6 Securing the sensor mat

The sensor mat is fastened with the lateral Velcro fasteners. These are joined together under the bed's patient surface. wissner-bosserhoff recommends that you fasten it to the seat section. Alternatively, you can also put the sensor mat on the lower third of the backrest. With the Easy-Clean patient surface, attachment is only possible on the seat section.

Recommended fastening:

Place the mat on the seat section. Fasten with Velcro fasteners.

Alternative fastening:

Place it on the lower third of the backrest. Fasten with Velcro fasteners.

Fasten any excess length of the sensor mat cable below the seat section with a cable tie.

Caution!

The sensor mats must not be bent, as this may stop them from functioning. For this reason, they may not be positioned in the transition zone between the seat section and backrest, since this would inevitably cause them to fold when the backrest is raised.

6.7 Connecting to the nurse call system or other signal receivers

All components required for connection (with the exception of any Y adapter that may be necessary) are supplied with the system. In order to ensure that the system functions smoothly with your nurse call system, certain things must be noted:

The adapter cable must be designed individually for the nurse call system installed in the house and must be ordered separately. Connect the system to the nurse call system by attaching the nurse call cable to the adapter cable through the XLR connector and the other end to the call system.

To allow you to continue using the handheld switch on your nurse call system and also to connect SafeSense®, two sockets are necessary. If your nurse call system does not have a socket, a Y-adapter is required. You can obtain one from your call system manufacturer. Please also ask the manufacturer about any risks and other information that should be noted.

If there is a Y cable for the corresponding call system, both connectors on the Y cable may need to be occupied (e.g. by a handheld switch).

wissner-bosserhoff is not responsible for the incorrect use of Y adapters or any possible damage to the call system that may result, and accepts no liability for such damage.

Attention!

wissner-bosserhoff accepts no liability for untransmitted or incorrectly transmitted signals to the nurse call system or other signal receivers. Before installing the SafeSense® Bed-Exit System, therefore, please contact the manufacturer to check the compatibility of your nurse call system and any other information that may need to be taken into account.

Calibrating a wireless receiver

The first time the wireless SafeSense® is connected to your nurse call system, the wireless receiver must be calibrated. This is done in three steps:


1. Press the button on the front of the wireless receiver. The light will start to blink.
2. Trigger an alarm by setting the time interval to 0 seconds and putting pressure on the sensor mat with your hand. After three seconds, remove your hand from the sensor mat. The blinking of the light on the wireless receiver will change.
3. Press the button on the front of the wireless receiver again. This is now calibrated to the SafeSense® system and will transmit the signal to the nurse call system.

Please note: The wireless receiver can be programmed to up to 30 senders. Therefore, for safety reasons, only one SafeSense® system may be calibrated to one wireless receiver at any time. If you wish to move the system to another location (e.g. another room), you will need to take the wireless receiver as well.

7. Operation of bed exit box

7.1 Switching on/off

Use the push-turn dial on the right side of the bed exit box to switch it on or off.

To switch the box on: press the push-turn dial for about 1 second

To switch the box off: press the push-turn dial for more than 3 seconds until the display turns off

7.2 Setting a time interval

Use the push-turn dial on the right side of the box to set the nurse call trigger time interval between 1 and 30 minutes. The tolerance between the time set and the actual time the call is triggered can be about 1 %.

➡ No time can be set when the timer is running (i.e. when the resident is not in bed). The timer can only be set if there is weight on the mat.

7.3 Display

The display on the top of the bed exit box shows the currently set delay for triggering the nurse call.

➡ If the bed exit box is switched off, the display shows the last time set when it is switched on again.

7.4 Underbed light

The Bed-Exit Box is equipped with an under-bed light, which supports the resident/patient at the nightly orientation. It gets activated, as soon as the resident/patient leaves the bed, regardless of whether a time interval is set, or not. As soon as the resident/patient returns to bed, the under-bed light switches off.

The under-bed light can also be controlled using the bed handset (from software version 1.5) of **nursing home beds**. However, this will only work if the Bed-Exit Box is activated. If the box is switched off, the under-bed light cannot be controlled with the handset. The underbed light is controlled with the handset as follows:

Activate the backrest with the following button:

Press the two arrow buttons (up and down) at the same time:

In this way, the under-bed light can be switched on or off.

- ➔ The underbed light cannot be switched off using the handset if it is on because of a call signal. In this case, it will only turn off if the Bed-Exit box is completely switched off or the resident returns to bed.
- ➔ Manual shifting operation with the hand control are only valid once. Once the resident leaves the bed or lying down, the automatic mode is active again.

8. Permissible mattresses and minimum resident weight

The sensor mat reacts to pressure relief, so observe the specified minimum weight of residents/patients in order for the system to reliably recognize bed egress. This depends on the type of mattress, the patient surface or the bed. The following charts show the weights above which the system reacts reliably in combination with the different mattresses and patient surfaces. Nevertheless, test the functionality with each individual resident/patient before use.

Permissible mattresses			Minimum resident weight
Mattress type	Mattress thickness	Articlenumber	for all types of patient surfaces
Universal mattress		50-0084	40 kg
ViskoMatt 2	14 cm	59-0143	40 kg
ViskoMatt 3	14 cm	50-0577	40 kg
ViskoMatt Komfort	14 cm	50-0317	40 kg
PrimaCare 10	14 cm	50-0705/ 50-0840	40 kg
PrimaCare 20	14 cm	50-0706	40 kg
CliniCare 10	14 cm	50-0707	40 kg
CliniCare 20	14 cm	50-0708	40 kg
CliniCare 30	14 cm	50-0709	40 kg
MicroMatt 7	14 cm	50-0565	40 kg

Attention!

If the resident is lighter than the specified weight, the sensor mat may not detect him/her leaving the bed, meaning that no signal will be sent to the nurse call system or any other signal recipient.

- ➔ Please consider that the sensor mats should not be used in combination with alternating pressure mattresses. Otherwise, the system cannot transmit reliable signals.
- ➔ The sensor mat is designed to be used with 90 cm wide patient surfaces. If the mat is used on a wider patient surface, take into account that some areas will not be covered by the mat. This might trigger the bed exit alarm even if the resident didn't leave the bed.

The data explicitly relates to wissner-bosserhoff mattresses. We offer no guarantee for these functions if other mattresses are used. If you are using a different mattress, please check its functionality and compatibility with the SafeSense® Bed-Exit System in advance.

9. Cleaning the sensor mat

The sensor mat, including the connection cable, must be cleaned regularly. This is ideally done by wiping it down with a disinfectant. The following cleaning agents can be used:

- Isopropanol (required acc. to EN 60601-1 art. 6.1)
- Denatured alcohol (required acc. to EN 60601-1 art. 6.1)
- Terralin
- Perform
- Microbac Forte
- Dismozon Pur
- Mikrozid
- Termosept

Caution!

The sensor mat has to be completely dry before it is used again. Do not clean the mat in the laundry or with water spray.

10. Maintenance and service

The SafeSense® Bed-Exit System requires minimal servicing. This is because, during product development, attention was paid to ensure that servicing work was reduced to a minimum, keeping operating costs low.

In daily use, however, experience shows that products are also sometimes handled carelessly, and rough handling can also contribute to them aging more quickly and to the wear of certain components, on which the manufacturer can have no direct influence.

For this reason, routine maintenance checks should be performed by the operator – this should also be done to guarantee the availability of the system, which is in the operator's own interest. We recommend servicing the Bed-Exit System once a year, together with the service of the nursing home bed.

Warranty or guarantee claims can only be enforced if the product has been serviced and maintained regularly and in accordance with the instructions for use and the operation and safety instructions in these instructions for use.

Only trained operators and maintenance technicians can ensure proper handling and use. wissner-bosserhoff's technical customer service offers servicing and necessary training on all aspects of the product.

If a function check, an inspection, measurement or servicing reveals serious defects that cannot be rectified, the product must be suspended from further use.

Our customer service is available to you for any inquiries and to provide training and servicing.

11. Spare parts

Only wissner-bosserhoff GmbH original spare parts may be used. Our customer service, sales and technical customer service provide information concerning spare part deliveries, etc. (for the address, see the Contact section).

Please request spare parts lists, current price lists and servicing instructions with exploded diagrams as required, supplying the data from the system's identification plate or the or the appropriate article number, order number and date of delivery, to wissner-bosserhoff, Technical Customer Service Department.

wissner-bosserhoff GmbH
customer service

Phone +49 2377 784-456
Fax +49 2377 784-150

Designation	Part No.
Universal Bed Exit box 2.1, cable	10-0812-0001
Universal Bed Exit box 2, wireless	10-0812-0002 (866 MHz, EU) 10-0813 (916 MHz, Australia, Canada)
Retaining plates	
Bed-Exit Box retaining plate for Duo 9	01-011193
Bed-Exit Box retaining plate for Q 7	01-011194
SafeLift retaining plate set, if SafeSense® also on the bed	02-001004
Sensor mat 2.1 80 x 20 cm	02-000876
Adapter cable	
Customer-specific adapter cable with on/off switch for connection to the nurse call system	04-1508-XXXX
Customer-specific adapter cable without on/off switch for connection to the nurse call system	04-1608-XXXX

Wireless receiver	
Call receiver (for connection to the NurseCall-system)	50-0681-0000 (EU) 50-0772-0000 (Australia) 50-0772-0001 (Canada)
Power plug call detector	01-010091
Plug-in receiver	50-0633
Miscellaneous	
Removable cable clip	01-003535

12. Disposal

This bed exit system falls within the scope of the EC Directive 2002/96/EC (WEEE). The operator must dispose of all components of the system appropriately in accordance with the provisions of the electrical and electronic equipment laws of your country.

When transferring the equipment to an industrial third party, you are bound by contract to point out that after the end of use, proper disposal must be undertaken or arranged. If you do not point this out to the third party, you are responsible for the proper disposal of the equipment after the end of use.

The metal and plastic parts that accumulate during service and repairs must be disposed of properly and professionally in accordance with the applicable laws and regulations. The electrical controls (Bed-Exit Box) in particular may only be disposed of through specialist firms or waste management facilities authorized to do this.

13. Warranty

Statutory warranty provisions apply.

This warranty covers all malfunctions and defects associated with materials and manufacture. Malfunctions and defects arising from incorrect handling and external influences are excluded from this. If there are grounds for justified complaints during the warranty period, these shall be dealt with free of charge. This warranty can be enforced with the sales receipt, which bears the date of purchase. Our general terms and conditions of delivery apply.

14. Technical data

Input voltage	24 – 30 V DC
Type of protection	IPX4
Protection class	2
Relative humidity	30% – 75%
Atmospheric pressure	800 hPa – 1060 hPa
Environment temperature	+10 °C – +40 °C
Manufacturer	 wissner-bosserhoff GmbH Hauptstraße 4 – 6 58739 Wickede (Ruhr) GERMANY Tel. +49 2377 784-0

15. Classification

According to Annex VIII, #13 of the Medical Device Directive (MDD) 2017/745, SafeSense® is a Class I medical device, as an accessory for medical beds.

Designation	Comment
MDR 2017/745	EU regulation for medical devices
Medical Device Law Implementation Act	MPDG (national implementation)
Directive 2014/53/EU	RED (Radio Equipment Directive)
DIN EN ISO 14971	Risk management applied to medical products
DIN EN 60601-1 (relevant sections)	Medical electrical appliances
DIN EN 60601-2-52 (applicable sections)	Medical beds
BfArM [German Federal Institute for Drugs and Medical Devices] recommendations	Recommendations of German Federal Institute for Drugs and Medical Devices
DIN EN 60529; VDE 0470-1	Type of protection by housing IP code (protection from humidity)
DIN VDE 0834	Standard for nurse call systems in hospitals and nursing homes

16. Electromagnetic compatibility

Table 1: Guidelines and manufacturer's declaration - electromagnetic emissions - for all ME equipment and ME systems

SafeSense® is designed to operate in an electromagnetic environment as specified below. The customer or user of SafeSense® should ensure that it is operated in such an environment.		
Interference measurement	Compliance	EMC guideline
RF emissions in accordance with CISPR 11	Group 1	SafeSense® uses RF energy exclusively for its internal function. Therefore its HF emission is very low and it is unlikely that neighbouring electronic equipment will be disturbed.
RF emissions in accordance with CISPR 11	Class [B]	SafeSense® is suitable for use in all types of institutions, including residential use and similar, directly connected to the public mains grid that also supplies buildings used for residential purposes.
Harmonics emission IEC 61000-3-2	[Complies]	
Voltage fluctuations/flicker emissions IEC 61000-3-3	[Complies]	

Table 2: Guidelines and manufacturer's declaration - electromagnetic immunity - for all ME equipment and ME systems

SafeSense® is designed to operate in an electromagnetic environment as specified below.			
Interference immunity checks	IEC 60601 test level	Compliance level	EMC guidelines
Electrostatic discharge (ESD) (IEC 61000-4-2)	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient-/burst (IEC 61000-4-4)	± 2 kV for power cables ± 1 kV for input and output lines	± 2 kV for power cables ± 1 kV for input and output lines	Grid power quality should be that of a typical commercial or hospital environment
Surge (IEC 61000-4-5)	± 1 kV voltage line-line ± 2kV voltage line-earth	± 1 kV voltage line-line ± 2 kV voltage line-earth	Grid power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines (IEC 61000-4-11)	< 5% U_T (> 95% dip in U_T) for ½ cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (95% dip in U_T) for 5 sec.	< 5% U_T (> 95% dip in U_T) for ½ cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (>95% dip in U_T) for 5 sec	Grid power quality should be that of a typical commercial or hospital environment. If the user requires SafeSense® to continue to function even in the event of interruptions to the power supply, we recommend connecting the Bed-Exit System to an uninterruptible power supply or a battery.
Magnetic field power frequency (50/60 Hz) (IEC 61000-4-8)	30 A/m	30 A/m	The magnetic field power frequency should be at levels characteristic for a commercial or hospital environment.
NOTE: U_T is the AC supply voltage prior to application of the test level			

Table 3: Guidelines and manufacturer's declaration - electromagnetic immunity - for all ME equipment and ME systems that are not life-sustaining

SafeSense® is designed to operate in an electromagnetic environment as specified below. The customer or user of SafeSense® should ensure that it is operated in such an environment.			
Interference immunity checks	IEC 60601-test level	Compliance level	EMC guidelines
Conducted RF (IEC 61000-4-6)	3 Vrms 150kHz bis 80MHz	3 Vrms	Portable and mobile communication devices should not be used when any closer to SafeSense®, including the cables, than the recommended separation distance calculated according to the applicable equation for the transmission frequency. Recommended separation distance $d=0,35*\sqrt{P}$ $d=0,35*\sqrt{P}$ 80 MHz bis 800 MHz $d=0,35*\sqrt{P}$ 800 MHz bis 2,7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol. 
	6 Vrms 150kHz bis 80MHz	6 Vrms	
Radiated RF (IEC 61000-4-3)	3 V/m 80MHz bis 2,7GHz	3 V/m	
	10 V/m 80MHz bis 2,7GHz	10 V/m	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			


- a Field strength from fixed transmitters such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot theoretically be predicted with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the field strength measured at the location where the SafeSense® Bed-Exit System is used exceeds the above compliance levels, the product should be monitored to verify that it is functioning as intended. If unusual performance is observed, additional steps may be required such as a change in orientation or a different location for the system.
- b For frequencies between 150 kHz and 80 MHz, the field strength should be lower than 3 V/m.

Table 4: Recommended safety distances between portable and mobile RF telecommunications equipment and ME equipment or ME system - for ME equipment or ME systems that are not life-supporting

Recommended separation distances between portable and mobile RF telecommunications devices and the SafeSense® Bed-Exit System.			
The SafeSense® Bed-Exit System is intended for operation in an electromagnetic environment where RF interference is controlled. The customer or SafeSense® user can therefore help to avoid electromagnetic interference by maintaining the minimum distance between portable and mobile RF telecommunications devices (senders) and the Bed-Exit System, depending on the communications device's power output as specified below.			
Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d=0,35*\sqrt{P}$	80 MHz to 800 MHz $d=0,35*\sqrt{P}$	800 MHz to 2,7 GHz $d=0,35*\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power of the transmitter in watts (W) according to the manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

17. Nameplate

Serial-No.: 0000000000-XXX

0000000000 = Order number

XXX = Running number per order

Article-No.: AA-AAAA VVVV

AA-AAAA = Article-No.

VVVV = Variant

PCB-Ver.: P.P

P.P = Actual board

Software-Ver.: S.S

S.S = Software-Version

18. Transport and storage

The same ambient conditions described in the technical data overview provided under section 14 apply when transporting and storing the bed exit system.

Especially during transport and storage, ensure that the sensor mat is not bent and that no heavy goods are on top of it.

19. Troubleshooting

	Errors/faults	possible cause	Solution
1.	No call is triggered and the under-bed light does not switch on	Bed-Exit box is switched off	Switch the box on
		Sensor mat cable is not secure	Plug the cable securely into the Bed-Exit Box
		Bed-Exit box is defective	Contact WiBo service
		No mains supply (bed with and without battery); battery is empty (bed with battery)	Plug the network cable into the socket
		Sensor mat is worn/defective	Replace sensor mat
2.	No call is triggered, but the under-bed light is on	The plug to the nurse call system is not secure	Check the plug connections, otherwise call WiBo service
		Time interval is not set correctly	Check and correct interval setting
		On/off switch on the adapter is off	Set the switch to "on"
3.	Call is triggered in the wrong room/incorrect calls/no calls are triggered	For wireless version: the sender is calibrated to the wrong receiver	Delete programming and recalibrate
4.	The call is triggered, but the under-bed light does not switch on	Light source defective	Contact WiBo service
		Bed-Exit box is defective	Contact WiBo service
5.	Calls are triggered incorrectly	The patient is not positioned correctly in the bed	Position patients correctly
		The sensor mat is not positioned correctly	Correct the position of the sensor mat and secure it in place
		The time interval set is too short	Adjust the time interval using the control dial
		Unsuitable mattress (too firm, too heavy)	Choose permissible mattress in accordance with the IFU

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meets the provisions of the medical device directive (MDR) 2017/745 and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device. /

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
when used as intended, complies with the requirements of the RED Directive 2014/53/EU. /

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complies with the essential requirements of the directive RoHS 2011/65/EU /

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Reproduction, also excerpts, only with prior permission of the publisher. All rights reserved.

Subject to changes due to technical developments.

All technical data are rated data and are subject to construction and manufacturing tolerances.

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