

Staebel Saddle Stool

Product Code: 111311

A height-adjustable, saddle stool featuring a comfortable padded seat.





This ergonomically designed saddle stool helps to reduce the discomfort of sitting for long periods and allows the user free range of upper body movement.

Product Features:

- High quality gas stool features a deluxe comfortable padded seat
- Five castor base offers greater stability and easy manoeuvrability
- Anti-bacterial upholstery

Specifications:

Product Code	111311
Height	560 - 630 mm
Seat Diameter	350 mm
SWL	150 kg