

Comflex Care Bahama Recliner

Product Code: 5022

Carefully designed to keep your patients or residents social at meal times, fitting in closely to many dining tables.



Our Complex Care Chair is a popular choice for many healthcare environments. It's available in gas or electric formats, easy to maneuver and simple to clean for infection control. It's been carefully designed to keep your patients or residents social at meal times, fitting in closely to many dining tables.

Product Features:

- Gas Assisted design
 - ‘Tilt in space’ seat and leg-rest action
 - Euro style reclining backrest with built-in Bahama style head wings
 - Elevating standard leg-rest
 - Incorporates pressure management Multicush gel seating with air cushion backrest system
 - Fixed lounge style armrest
 - Includes push handle
 - 180kg Safe Working Load
 - Multi stretch material
 - Manufactured in Australia
- *Available in a huge range of upholstery colours



Specifications:

Product Code	5022
Seat Width	550 mm
Seat Depth	520 mm
Seat Height	550 mm
Overall Width	750 mm
Overall Depth (Upright/Flat)	830 / 1800 mm
Overall Height	1200 mm
SWL	180 kg

