

Regency Result3 Electric Recliner

Product Code: 5025

Packed full of quality features, this sleek recliner boasts a modern design appealing to wards, aged care facilities and clinics.



Packed full of features and with a 200kg weight capacity, this sleek recliner boasts a modern design appealing to wards, aged care facilities and clinics. The fixed lounge-style armrest has been uniquely designed to fit right up to most dining tables to increase social interaction. Created with infection control and easy-clean in mind, this chair is a great addition to a range of healthcare settings.

Product Features:

- Electric ‘tilt in space’ seat and legrest action
- Electric Euro style reclining backrest & elevating legrest
- Incorporating pressure management Multicush gel seating and air cushion backrest system
- Auto lay flat function and auto chair function
- Including push handle
- 4x 175mm swivel braking castors
- Inclusive seams for less pinch points
- Sleek, modern with all vinyl design and fixed lounge-style armrest
- Manufactured in Australia



Specifications:

Product Code	5025
Seat Width	550 mm
Seat Depth	520 mm
Seat Height	550 mm
Backrest Height	730 mm
Overall Width	810 mm
Overall Height	1280 mm
SWL	200 kg

