

Regency R2700 Flotation Chair

Product Code: 5030

This Regency R2700 Series Care Chair is designed to be carer and client friendly.





This Regency R2700 Series Care Chair is designed to be carer and client friendly by incorporating the original tilt seat system into a multi-positioning, multi-functional and manoeuvrable care chair. This particular model has a narrow seat width and shorter seat depth for the small adult or paediatric client.

Product Benefits:

- Superior range of tilt and recline in comparison with many products on the market
- Low seat height
- Pressure management seating – gel cushion in seat, air in backrest
- Forward seat tilt for ease to stand, legrest angle under seat
- Elevating legrest complete with fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Drop down armrests and swingaway headwings for access
- Adjustable head rest cushion
- 125mm castors, braking on rear, directional lock on front
- Quality Australian made construction
- Various colours available for arm rests and side rests

Specifications:

Product Code	5030
Seat Width	410 mm
Seat Depth	420 mm
Seat Height (Tilt/Normal)	520 mm / 550 mm
Backrest Height	700 mm
Armrest Height	175 mm
Legrest Length	470 – 550 mm
Overall Width	660 mm
Overall Depth	1050 mm
Overall Height	1280 mm
SWL	180 kg