


Commode Transfer Bench

Product Code: 1019

The Commode Transfer Bench has all the features of the standard transfer bench with the extra benefit of padded seat.



Product Features:

- Two piece padded seat allows water to drain away from seating area
- Ideal solution for safe bathroom transfers
- Innovative designs for shower, bath and over toilet
- Bridges many difficulties with hip baths as well as fibreglass enclosures
- Tool free assembly - snap on/off for left/right hand positioning
- Excellent for amputees or patients who must first sit outside the bath before placing legs into the bath
- Side arms aids in lateral transfer & positioning
- Easily adjustable backrest & legs
- Rack that allows a bowl to slide under for commode use

Specifications:

Product Code	1019
Seat Width	630 mm
Seat Depth	405 mm
Seat Height	460-560 mm
Bath Rim Clear	340-440 mm
Seat-Back Height	380 mm
Overall Width	730 mm
Overall Depth	600 mm
SWL	125 kg