

Staebel Transfer Bench

Product Code: 10TB

Designed for individuals with restricted mobility, the Standard Transfer Bench provides safety & independence for all users.



1018J



1018J-1



Designed for individuals with restricted mobility, the Standard Transfer Bench provides safety and independence while transferring in and out of the bath tub. This durable aluminium transfer bench is adjustable in height.

Product Features:

- 3 piece plastic non-slip seat with drain holes
- Hand rail for added support while transferring
- Reversible backrest to suit every situation
- 4 Screws to attach backrest to frame
- Welded frame
- Available in plastic or padded design
- Safe Working Load of 140kg

Specifications:

Product Code	1018J	1018J - 1
Seat Width	670 mm	670 mm
Seat Depth	410 mm	410 mm
Seat Height	445 - 570 mm	460 - 585 mm
Bath Rim Clear	370 - 495 mm	390 - 490 mm
Seat-Back Height	410 mm	410 mm
Overall Width	790 mm	790 mm
Overall Depth	540 mm	540 mm
SWL	140 kg	140 kg