

Staebel Shower Trolley

Product Code: 1209ST

The Staebel Shower Trolley features a one-section top and electrically adjustable positions to help with showering.



Product Features:

- Shower trolley with one-section top
- Central braking system
- Easy steering system (for ease of manoeuvring trolley's front left castor can be locked in direction of steering)
- Foldable side rails with a safety release mechanism (material – powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 270 mm)
- Mattress integrated drain with 645 – 2000 mm long flexible hose
- Water-resistant head pillow – semi – cylinder
- Safe Working Load 170kg

Specifications:

Product Code	1209ST
Shower trolley top Width	650 mm
Shower trolley top Length	1900 mm
Shower trolley top Height	540 - 970 mm
Height with Castors	600 - 1030 mm
Trolley top Inclination	0° - 5°
Length	2070 mm
Width	1070 mm
SWL	170 kg