

Staelbel Shower Trolley with Backrest

Product Code: 1210ST

The Staelbel Shower Trolley with backrest features a two-section top and electrically adjustable positions to help with showering.



Product Features:

- Shower trolley with a two-piece trolley top
- Adjustable 5° trolley top steered via gas springs
- Foldable side rails with safety mechanismHydraulically or electrically adjustable height of the trolley top in the range between 54cm -97cm Central braking system
- Easy steering system (for ease of manoeuvring the trolley's front left castor can be locked in the direction of steering)
- Metal sheet base cover
- Mattress with flexible drain hose, length from 70cm to 200cm
- Water-resistant head pillow semi-cylinder
- Metal frame powder-coated with antibacterial paint

Specifications:

Product Code	1210ST
Shower trolley top Width	810 mm
Shower trolley top Length	1900 mm
Trolley top Inclination	0° - 5°
Head Section Inclination	25°
Electric Hi-Lo	540 - 970 mm
Length	2090 mm
Width	810 mm
SWL	200 kg

Issue Date: 150324
Revision Date: 190525 | V2