

Staebel Electric Shower Trolley with Backrest

Product Code: 1211ST

The Staebel Shower Trolley – Electric with Backrest features a versatile two-section top and electrically adjustable positions to help with showering in healthcare environments.



Product Features:

- Adjustable 5° trolley top inclination
- Electrically lowering/raising of trolley top via console
- Central braking system
- Easy steering system (for ease of manoeuvring trolley's front left castor can be locked in direction of steering)
- Foldable side rails with a safety release mechanism (material – powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 27 cm)
- Mattress integrated drain with 645 – 2000 mm long flexible hose
- Water-resistant head pillow

Specifications:

Product Code	1211ST
Shower Trolley Top Width	650 mm
Shower Trolley Top Length	1900 mm
Shower Trolley Top Height	540 – 970 mm
Trolley Head Section Inclination	45°
Width	810 mm
Length	2070 mm
Castors	150 mm
SWL	170 kg