

Staelbel Electric Profiling Shower Trolley

Product Code: 1213ST

The Staelbel Profiling Shower Trolley has been designed to facilitate shower bathing and hygiene routines, for users and carers with a wide range of postural and cognitive challenges.



Product Features:

- Shower trolley with four-section top
- Electrically adjustable with Trendelenburg and, Anti-Trendelenburg functions
- Foldable side rails with a safety release mechanism (material - powder coated stainless steel)
- 150 mm with central braking castors
- Easy steering system (for ease of manoeuvring trolley's front left castor can be locked in direction of steering)
- Metal base cover
- Water-resistant shower mattress
- Mattress integrated drain with 100cm long flexible hose
- Water-resistant head pillow – semi-cylinder

Specifications:

Product Code	1213ST
Shower Trolley Top Width	650 mm
Shower Trolley Top Length	1900 mm
Trolley Head Section Inclination	0° - 60°
Electric Hi-Lo	520 - 920 mm
Length	2070 mm
Width	810 mm
SWL	200 kg

Issue Date: 160124
Revision Date: 190525 | V2