

Staebel Junior Shower Trolley

Product Code: 1214ST

The Staebel Junior Shower Trolley features a one-section top and electrically adjustable positions to help with showering, specifically for paediatric or smaller users.



Product Features:

- Adjustable 5° trolley top inclination
- Hydraulic lowering/raising of trolley top via a console from 54 cm to 97 cm
- Central braking system for safety
- Foldable side rails with a safety release mechanism (material – powder coated stainless steel)
- Easy steering system (for ease of manoeuvring trolley's front left castor can be locked in direction of steering)
- Metal base cover
- Water-resistant shower mattress (height: 27 cm)
- Mattress integrated drain from 645 – 2000 mm long flexible hose
- Water-resistant head pillow – semi-cylinder

Specifications:

Product Code	1214ST
Shower Trolley Top Width	600 mm
Shower Trolley Top Length	1900 mm
Length	1670 mm
Shower Trolley Top Height	540 - 970 mm
Width	760 mm
Castors	150 mm
SWL	170 kg

Issue Date: 150324
Revision Date: 281025 | V2