

Staebel Bariatric Shower Trolley

Product Code: 1216ST

The Staebel Bariatric Shower Trolley features a one-section top and electrically adjustable positions to help with showering, specifically for bariatric users.



Product Features:

- Hydraulic lowering/raising of the stretcher top via a console lift
- Castors with a diameter of 150 mm
- Central braking easy steering system (for ease of manoeuvring the trolley's front left castor can be locked in the direction of steering directional wheel)
- Foldable side rails with a safety release mechanism (material – powder coated stainless steel)
- Metal base cover
- Water resistant shower mattress (height: 270mm)
- Mattress integrated drain with a 645- 2000mm long flexible hose
- Water resistant head pillow – semi – cylinder

Specifications:

Product Code	1216ST
Shower Trolley Top Width	900 mm
Shower Trolley Top Length	2000 mm
Shower Trolley Top Height	540 - 970 mm
Height with Castors	600 - 1030 mm
Trolley Top Inclination	0° - 5°
Length	2170 mm
Width	1070 mm
SWL	250 kg