

# Staebel Folding Shower Chair

Product Code: 1350

The Staebel Folding Shower Chair, with the ability to fold the frame, is the solution to those needing a Shower Chair but live in multiple locations or travel.



## Product Features:

- Foldable Design for easy travel or storage
- Non-slip textured surface
- Moulded Plastic Seat
- Front cut-out for personal hygiene
- Adjustable seat height to suit individual needs
- Rust-proof epoxy-coated aluminium frame
- Safe Working Load 100kg

## Specifications:

|                    |              |
|--------------------|--------------|
| Product Code       | 1350         |
| Seat Width         | 410 mm       |
| Seat Depth         | 400 mm       |
| Seat Height        | 455 - 520 mm |
| Width Between Arms | 470 mm       |
| Overall Width      | 530 mm       |
| Overall Depth      | 455 mm       |
| SWL                | 100 kg       |