

Staebel Shower Stool with Arms & Back

Product Code: 1569

This Shower Stool with Arms and Backrest, provides security when entering or exiting the shower and assist the user on and off the seat.



This Staebel Shower Stool with Arms and Backrest, provides security when entering or exiting the shower and assist the user on and off the seat.

Product Features:

- Integrated arms assist transfers
- Non-slip rubber tips
- Electrophoretic surface treatment to prevent rusting
- Powder coated steel frame
- Compact size for the smaller shower recess
- Aluminium construction for long life
- Padded seat and backrest for user comfort
- Adjustable seat height

Specifications:

Product Code	1569
Seat Width	370 mm
Seat Depth	290 mm
Seat Height	520 - 670 mm
Width Between Arms	460 mm
Armrest Height (From Floor)	640 - 840 mm
Overall Width	550 mm
Overall Depth	450 mm
SWL	125 kg