

Bariatric Shower Chair with Padded Seat

Product Code: 1579

The Bariatric Shower Chair is designed to assist in user showering with a comfortable and durable design.



Product Features:

- Padded seat and curved backrest provide comfort and postural support
- Powder coated steel frame
- Zinc treated to prevent rusting
- Made in Australia
- 200kg SWL
- Suitable for Bariatric Users

Specifications:

Product Code	1579
Seat Width	600 mm
Seat Height	520 mm
Seat Depth	450 mm
SWL	200 kg