

Aluminium 3-in-1 Over Toilet Frame

Product Code: 1752

Designed with versatility and user comfort in mind, this frame serves as a multi-functional aid to support individuals with limited mobility.



Product Features:

- The 3-in-1 frame can be used as an over-toilet aid, a standalone commode, or a shower chair, making it a versatile solution for various personal care needs.
- Crafted from high-quality aluminium, the frame is lightweight yet sturdy. Its rust-resistant properties ensure long-lasting durability, even in humid bathroom environments.
- Adjustable Height: Easily adjust the frame's height to accommodate different users and ensure maximum comfort and safety.
- Comes complete with a splashguard and a removable bowl with a lid, making it ideal for use as a commode. The bowl can be removed easily and cleaned.

Specifications:

Product Code	1755
Width Between Arms	450 mm
Seat Height	475 - 625 mm
Overall Width	560 mm
Armrest Height	700 - 870 mm
Footprint (W x D)	560 x 590 mm
SWL	125 kg