

Staebel Over Toilet Frame

Product Code: 17OTF

Designed to assist the user who has difficulty sitting down or standing up by adjusting the toilet seat height.



Product Features:

- Zinc treated, epoxy coated frame with armrests and plastic seat
- Removable splash guard included
- Lift up seat that removes for cleaning
- Ideal for those with hip or knee injuries or who require an increased seat height
- Unique stacking design saves space when storing more than one frame
- Multi D leg adjustment ensures stability on uneven surfaces
- Push button height adjustment
- Non slip guides on base of legs
- Safe Working Load of 150kg

Specifications:

Product Code	1750J
Width Between Arms	500 mm
Seat Width	370 mm
Overall Width	590 mm
Seat Depth	430 mm
Seat Height	470 - 610 mm
Footprint (WxD)	590 x 620 mm
Armrest Height	750-855 mm
SWL	150 kg

Issue Date: 010222
Revision Date: 190226 | V4