

# Toilet Surround

Product Code: 190

The Toilet Surround fits around the toilet and provides assistance in getting on and off the toilet.



## Product Features:

- Heavy Duty Design
- Ergonomic shaped side frame makes it easy to stand
- Non- slip feet
- Push button for height adjustment
- Steel manufactured zinc treated and epoxy coated making it both strong and durable

## Specifications:

Product Code	1900	1901	1902
Width Between Arms	460 mm	540 mm	620 mm
Arm Height	685 mm	685 mm	685 mm
Overall Width	540 mm	605 mm	670 mm
Overall Depth	490 mm	520 mm	520 mm
SWL	125 kg	130 kg	150 kg