

Juvo Paediatric Attendant Propelled Stainless Steel Shower Commode

Product Code: 190JAS

Designed to assist in transportation of the user for showering and toileting, specifically paediatric or smaller users.




This Juvo Paediatric Shower Commode assists transportation of the user for showering and toileting, specifically paediatric or smaller users. It is made from durable total stainless steel frame and can be positioned over the toilet.

The patented TearSafe design ensures safety by eliminating protrusions, and the 125mm stainless steel castors, equipped with four brakes, enhance stability. Additional features include swing-back, independently adjustable armrests for easier side transfers and patient access, along with two-position height-adjustable armrests for added comfort.

Product Features:

- Adjustable Backrest design – adjust the backrest angle from 90 to 120 degrees to improve client comfort and access for the carer
- Depth adjustable seat – adjust the seat depth to suit smaller client by up to 75mm
- Height adjustable seat up to 90mm
- 125mm stainless steel braking castors
- Self-propelled option available – frame ready to accept self-propelling kit
- Tested and certified to Australian Standards

Please note that the seat for this product is sold separately (See colours attached)



Lemon Twist

Lime

Lipstick

Regal

Emerald

Specifications:

Product Code	1903J	1905JDNU
Seat Width	360 mm	410 mm
Width Between Arms	390 mm	440 mm
Anterior Seat Height	595 mm	595 mm
Posterior Seat Height	565 mm	565 mm
Overall Width	440 mm	490 mm
SWL	136 kg	136 kg