

Juvo Self Propelled Shower Commode with Swing Away Footrests

Product Code: 1941J

The ultimate in strength and quality, it is designed to assist in transportation of the user for showering and toileting





The Juvo Self Propelled Mobile Shower Commode is the ultimate shower commode to assist in the transportation of the user for showering and toileting.

Product Features:

- Self propelled
- Stainless steel frame for longevity
- Flip up armrests for ease of side transfers
- Slatted backrest design which facilitates easy cleaning
- Safety armrests can be flipped so that they are positioned in front of the user to eliminate slipping forward
- Swingaway adjustable footplates with unique design that reduces skin tear risk
- Tested and certified to Australian Standards

Please note that the seat for this product is sold separately (See colours attached)

Specifications:

Product Code	1941J
Seat Width	410 mm
Width Between Arms	440 mm
Overall Width	610 mm
Anterior Seat Height	595 mm
Posterior Seat Height	565 mm
Armrest Height	200 - 240 mm
SWL	200 kg

Lemon Twist

Lime

Lipstick

Regal

Emerald