

# RAZ CAT Compact Mobile Shower Commode

Product Code: 1946

The RAZ-CAT Mobile Shower Commode Chair has all of the features of the RAZ-AT but in a more compact size.



The RAZ-CAT Mobile Shower Commode Chair has all of the features of the RAZ-AT but in a more compact size. For pediatrics, adolescents or smaller adults, the RAZ-CAT provides a fit that is better suited to physically smaller clients. The depth-adjustable back provides growth capability, an essential component for pediatrics. Designed and manufactured to the highest standards of the discerning healthcare professional.

## Product Features:

- Base and seat frames retain the 18" width of the Raz-AT so that the chair can still be positioned over a toilet.
- Base frame has been shorted by 3", which makes the chair more maneuverable.
- The standard seat has an elongated aperture that works with seat depths from 11" to 16".
- Seat frame depth has been reduced and a depth-adjustable back offers depths from 12" to 15" (to front of frame) or up to 17" with the front of the seat overhanging the seat frame.

## Specifications:

Product Code	1946-16	1946-18	1946-20
Frame Width	460 mm	460 mm	460 mm
Overall Width	<b>460 mm</b>	<b>550 mm</b>	<b>580 mm</b>
Overall Length	810 mm	810 mm	810 mm
Seat Width	460 mm	460 mm	460 mm
Seat Depth	280 – 400 mm	280 – 400 mm	280 – 400 mm
Armrest Height (From Seat)	140 – 200 mm	140 – 200 mm	140 – 200 mm
Seat to Floor Height	550 – 650 mm	550 – 650 mm	550 – 650 mm
Product Weight	25 kg	25 kg	25 kg
SWL	115 kg	115 kg	115 kg