

Raz AP Adjustable Height Shower Commode

Product Code: 1947

Available with so many adjustments and options available, the Raz-AP is a true prescriber-friendly, mobile shower chair.



The Raz-AP mobile shower commode chair comes standard with 4" of height adjustment, an adjustable-tension fabric backrest, 5" dual-locking casters, flip-up padded armrests, and a commode pan. The easy-to-clean, polished stainless steel frame is designed for easy side and front access by the user or attendant. Only Raz offers the Ischial & Pelvic Alignment System (IPAS), which allows the commode seat aperture to be comfortably centered around the ischial tuberosities for clients with a posterior or anterior rotated pelvis.

Product Features:

- Adjustable seat height
- Padded arm rests
- Adjustable angle backrest with adjustable tension
- 4 x Dual-locking (swivel and roll) castors
- Ischial & Pelvic Alignment System (IPAS), which allows the commode seat aperture to be comfortably centered around the ischial tuberosities.
- Adjustable / Quick-Release Seat
- Swing away legrests
- Safe Working Load 158kg

Specifications:

Product Code	1947
Frame Width	460 mm
Overall Width	570 mm
Overall Length	860 mm
Seat Width	460 mm
Seat Depth	400 mm / 480 mm
Clearance Height (Under Seat) (with / without pan brackets)	390 - 500 mm 410 - 510 mm
Armrest Height (From Seat)	170 - 240 mm
Seat to Floor Height	510 - 610 mm
Product Weight	17.25 kg
SWL	158 kg