

RAZ SP Adjustable Height Shower Commode

Product Code: 1949

Available with so many adjustments and options available, the Raz-SP is a true prescriber-friendly, mobile shower chair.



The Raz-SP mobile shower commode chair comes standard with 4" of height adjustment, an adjustable-tension fabric backrest, 22" rear wheels, 5" dual-locking casters, flip-up padded armrests, and a commode pan. The easy-to-clean, polished stainless steel frame is designed for easy side and front access by the user or attendant.

Only Raz offers the Ischial & Pelvic Alignment System (IPAS), which allows the commode seat aperture to be comfortably centered around the ischial tuberosities for clients with a posterior or anterior rotated pelvis.

Product Features:

- Adjustable seat height
- Flip Up Padded arm rests
- Adjustable angle backrest with adjustable tension
- 4 x Dual-locking (swivel and roll) castors
- Ischial & Pelvic Alignment System (IPAS), which allows the commode seat aperture to be comfortably centered around the ischial tuberosities.
- Multi-Position Axle Plates
- Adjustable / Quick-Release Seat
- Locking gas struts to facilitate tilting

Specifications:

Product Code	1949-16	1949-18
Frame Width	460 mm	460 mm
Overall Width	670 mm	670 mm
Overall Length	990 mm	990 mm
Seat Width	460 mm	460 mm
Seat Depth	405 mm / 480 mm	405 mm / 480 mm
Armrest Height (From Seat)	180 - 250 mm	180 - 250 mm
Seat to Floor Height	510 - 610 mm	510 - 610 mm
Product Weight	21.3 kg	22 kg
SWL	160 kg	160 kg