

# RAZ Zum Folding AP Mobile Shower Commode

Product Code: 19569

The Züm-AP (Attendant Propel) is the first folding frame mobile shower commode chair that offer the rigidity that you would expect from a premium product.



Only Raz offers the Ischial & Pelvic Alignment System (IPAS), which allows the commode seat aperture to be comfortably centered around the ischial tuberosities for clients with a posterior or anterior rotated pelvis.

## Product Features:

- Crossbrace with low pivot point and forward arch allows Züm to be positioned fully over the toilet
- Crossbrace lock-down system and rear seat support tube provide unparalleled rigidity
- 18"W aluminum frame – fits over most standard toilets
- Unrestricted front and side access
- 17 standard seats + custom (with IPAS)
- Removable back support
- Flip-up arm supports with optional lock kit
- 2" fore/aft rear wheel adjustment (fine tunes center-of-mass position for SP model)
- Quick-Release Axles are standard

## Specifications:

<b>Product Code</b>	195769
<b>Frame Width</b>	470 mm
<b>Overall Width</b>	280 mm
<b>Overall Length</b>	870 mm
<b>Seat Width</b>	460 mm
<b>Clearance Height (Under Seat)</b> (with / without pan brackets)	430 / 450 mm
<b>Armrest Height (From Seat)</b>	200 mm
<b>Seat to Floor Height</b>	550 mm
<b>Product Weight</b>	20 kg
<b>SWL</b>	159 kg