

RAZ Zum Folding SP Mobile Shower Commode

Product Code: 19570

The Züm-SP (Self Propel) is the first folding frame mobile shower commode chair that offer the rigidity that you would expect from a premium product.



Only Raz offers the Ischial & Pelvic Alignment System (IPAS), which allows the commode seat aperture to be comfortably centered around the ischial tuberosities for clients with a posterior or anterior rotated pelvis.

Product Features:

- Crossbrace with low pivot point and forward arch allows Züm to be positioned fully over the toilet
- Crossbrace lock-down system and rear seat support tube provide unparalleled rigidity
- 18"W aluminum frame – fits over most standard toilets
- Unrestricted front and side access
- 17 standard seats + custom (with IPAS)
- Removable back support
- Flip-up arm supports with optional lock kit
- 22" or 24" wheels available for the Züm-SP model
- 2" fore/aft rear wheel adjustment (fine tunes center-of-mass position for SP model)
- Quick-Release Axles are standard

Specifications:

Product Code	19570
Frame Width	470 mm
Overall Width	640 mm
Overall Length (22" / 24" Wheels)	1020 / 1040 mm
Seat Width	460 mm
Clearance Height (Under Seat) (with / without pan brackets)	430 / 450 mm
Armrest Height (From Seat)	200 mm
Seat to Floor Height	550 mm
Product Weight	26 kg
SWL	159 kg