

# RAZ AT600 Heavy Duty Mobile Shower Commode

Product Code: 19572

The RAZ-AT600 – Heavy Duty Mobile Tilt Shower Commode Chair with stainless steel frame and components, 25° tilt range, four locking 6" (150mm) casters and commode pan.



Raz-AT600 – Heavy Duty Mobile Shower Commode Chair (up to 605 lbs. maximum) with stainless steel frame and components, 25° tilt range, four locking 6" (150mm) casters and commode pan.

The specially engineered, Weight Assist Spring System (WASS) is available in two weight ranges and is adjustable to accommodate the user's weight and Center of Mass (CoM) for tilt operation with minimal effort by the attendant. Headrest pad (10" x 4"), stainless steel headrest fixture, Tilt Assist Pedal and HD Footrests are included.

## Product Features:

- WASS – Weight Assist Spring System: this counter balances the user's weight to facilitate tilt. Unlike battery-powered MSCCs, the AT600 never needs charging and is always ready for use!
- 25° tilt range
- Adjustable back angle (90° / 97°)
- Back frame adjustable fore / aft (3" range)
- Gluteal shelf accommodation
- Height-adjustable
- 6" dual-locking casters
- Stainless steel frame



### Specifications:

Product Code	19572-22	19572-24	19572-26	19572-28	19572-30
Frame Width	560 mm	600 mm	660 mm	710 mm	760 mm
Overall Length	960 - 990 mm	960 - 990 mm	960 - 990 mm	960 - 990 mm	960 - 990 mm
Seat Width	571 mm	611 mm	671 mm	721 mm	771 mm
Seat Depth	460 mm	460 mm	460 mm	460 mm	460 mm
Armrest Height (From Seat)	240 mm	240 mm	240 mm	240 mm	240 mm
Seat to Floor Height	550 - 580 mm	550 - 580 mm	550 - 580 mm	550 - 580 mm	550 - 580 mm
Product Weight	37 kg / 38 kg	37 kg / 38 kg	37 kg / 38 kg	37 kg / 38 kg	37 kg / 38 kg
SWL	275 kg	275 kg	275 kg	275 kg	275 kg