

# RAZ AP Heavy Duty Mobile Shower Commode

Product Code: 1966

The Raz-AP600 will provide maximum comfort and maneuverability for users who require an attendant propelled heavy-duty shower commode chair.



The Raz-AP600 mobile shower commode chair features the Multi-Position Back (with center of gravity adjustment) and comes standard with 4" of height adjustment, an adjustable-tension fabric backrest, flip-up padded armrests, a commode pan, and a large molded soft foam seat. The electro-polished stainless steel frame is engineered for heavy-duty use and easy access for the user or attendant.

Only Raz offers the Ischial & Pelvic Alignment System (IPAS), which provides 2" of seat depth adjustment to accommodate a variety of body sizes and shapes.

## Product Features:

- Gluteal shelf accommodation
- Unrestricted front and side access
- Adjustable back angle (90° / 97°)
- Back frame adjustable fore / aft (2" range)
- 5" dual-locking casters (6" optional)
- Stainless steel frame
- SWL of 275kg

## Specifications:

Product Code	1966
Overall Width	670 mm
Overall Length	910 - 990 mm
Seat Width	580 mm
Seat Depth	460 mm
Armrest Height (From Seat)	240 mm
Seat to Floor Height	480 - 580 mm
Product Weight	25 kg / 6.7 kg
SWL	275 kg