

RAZ ART Tilt Mobile Shower Commode

Product Code: 1979

The Raz-ART (Attendant Rotational Tilt) mobile shower commode chair tilts up to 40° to provide extra pressure reduction and comfort.



The Raz-ART (Attendant Rotational Tilt) mobile shower commode chair tilts up to 40° to provide extra pressure reduction and comfort. The unique frame design utilizes a patented, dual-arc, orbital tilt system that allows for the smallest footprint and lowest tilt effort of any chair on the market. Patent-pending tilt-locking system of the Raz-ART eliminates cable “looping,” which is prevalent on other designs. The Raz-ART mobile shower commode chair comes standard with the Ischial & Pelvic Alignment System (IPAS), adjustable-tension back upholstery, molded headrest and fixture, 5” dual-locking casters, flip-up padded armrests, and a commode pan.

Product Features:

- Ischial & Pelvic Alignment System (IPAS)
- Adjustable-tension back upholstery
- Molded headrest and fixture
- 5” dual-locking casters
- Flip-up padded armrests, and a commode pan.
- The easy-to-clean, aluminum frame fits over most toilets and is designed for easy side and front access by the user or attendant.

Specifications:

Product Code	1979
Overall Width	570 mm
Overall Length (Upright / Tilted)	890 / 1280 mm
Seat Width	460 mm
Clearance Height (Under Seat) (with / without pan brackets)	380 - 560 mm 410 - 590 mm
Armrest Height (From Seat)	230 mm
Tilt Range	Tilted to 40°
Seat to Floor Height	510 mm - 585 mm / 610 - 685 mm
Product Weight	31 kg
SWL	160 kg