

Juvo Attendant Propelled Height Adjustable Shower Commode

Product Code: 19AAPSCA

The ultimate shower commode to assist in the transportation of the user for showering and toileting.



The Juvo Height Adjustable Attendant Propelled Commode is a versatile and durable solution designed to assist with patient transportation for showering and toileting. The commode can be conveniently positioned over a toilet and is built with a robust, total stainless steel frame for longevity. The patented TearSafe design ensures safety by eliminating protrusions, and the 125mm stainless steel castors, equipped with four brakes, enhance stability.

Additional features include swing-back, independently adjustable armrests for easier side transfers and patient access, along with two-position height-adjustable armrests for added comfort.

Product Features:

- Height adjustable up to 90mm front and rear axle
- Reinforced seat/back frame
- Swing away, removable and height adjustable legrests with flip-up footplates
- Patented design, no protrusions
- 125mm stainless steel braking castors
- Swing back independent armrests to assist in side transfer and patient access
- Safety arms provide user security & safety
- Stainless steel pan rack guides
- Tested and certified to Australian Standards

Please note that the seat for this product is sold separately (See colours attached)



Lemon Twist

Lime

Lipstick

Regal

Emerald

Specifications:

Product Code	1906J	1913J	1918J	1919J	1933J
Seat Width	410 mm	460 mm	510 mm	560 mm	660 mm
Anterior Seat Height	505 - 595 mm	505 - 595 mm	505 - 595 mm	505 - 595 mm	505 - 595 mm
Posterior Seat Height	475 - 565 mm	475 - 565 mm	475 - 565 mm	475 - 565 mm	475 - 565 mm
Width Between Arms	440 mm	490 mm	540 mm	590 mm	690 mm
Overall Width	490 mm	540 mm	590 mm	640 mm	740 mm
SWL	200 kg	200 kg	300 kg	300 kg	300 kg