

Juvo Attendant Propelled Shower Commode with Swingaway Legrests

Product Code: 19JASA

The ultimate in strength, quality and stability this Juvo mobile shower chair is built for heady duty use.



The ultimate in strength, quality and stability the Juvo Attendant Propelled Shower Commode with Swingaway Legrests is built for heavy use. This model is a mobile shower commode with swing away footrests. The long wheelbase stops forward tipping and has proven to be an effective solution in all areas from home to hospital.

Product Features:

- **Quality stainless steel frame for longevity in a demanding environment** - This product features a durable stainless steel frame that has been specifically designed to withstand the rigours of demanding environments. The use of high-quality stainless steel ensures resistance to wear, corrosion, and the general effects of daily use, making it an ideal choice for long-term reliability in care or medical settings.
- **Armrests can be flipped so that they are positioned in front of the occupant to reduce the risk of slipping forward** - The armrests on this product are designed to provide enhanced safety and adaptability. With a simple flipping mechanism, the armrests can be moved to the front of the occupant. This innovative feature helps minimise the risk of the user slipping forward, adding an extra layer of security during use.
- **Smooth, slatted backrest design to facilitate easy cleaning** - The backrest incorporates a smooth, slatted design that is both functional and hygienic. This design allows for effortless cleaning, ensuring that the product can be maintained to high cleanliness standards with minimal effort. Its practicality makes it particularly well-suited to environments where hygiene is a top priority.
- **Swing-away, removable, and height-adjustable legrests** - Designed for versatility and user comfort, the legrests can swing away for easy transfers, be completely removed when not needed, and are height adjustable to provide personalised support and ergonomic positioning.
- **This product is tested and certified to Australian Standards**
- **Stainless steel 125mm castors (4 braking)**

Please note that the seat for this product is sold separately
(See colours attached)



Lemon Twist

Lime

Lipstick

Regal

Emerald

Specifications:

Product Code	1905J	1915J	1916J	1917J	1935J
Seat Width	410 mm	460 mm	510 mm	560 mm	610 mm
Anterior Seat Height	595 mm	590 mm	590 mm	590 mm	590 mm
Posterior Seat Height	505 mm	565 mm	565 mm	565 mm	565 mm
Width Between Arms	440 mm	465 mm	520 mm	570 mm	660 mm
Overall Width	490 mm	540 mm	610 mm	640 mm	710 mm
SWL	200 kg	200 kg	300 kg	300 kg	300 kg