

Juvo Self Propelled Folding Mobile Shower Commode

Product Code: 19SPSCA

The ultimate shower commode to assist in the transportation of the user for showering and toileting.



Product Features:

- Assists transportation of the patient for showering and toileting
- Can be positioned over the toilet
- Durable total stainless steel frame
- Removable seat
- Reinforced seat/back frame
- Swing away, removable & height adjustable legrests with flip-up footplates.
- Patented design, no protrusions
- Now with braking stainless steel 125mm castors
- Swing back independent armrests to assist in side transfer and patient access
- Two position height adjustable armrests
- Safety arms provide user security & safety
- Stainless steel pan rack guides
- Max User Weight 200kg
- Australian Standards Approved
- Patented folding design into standard suitcase size for easy travel

Please note that the seat for this product is sold separately (See colours attached)



Lemon Twist

Lime

Lipstick

Regal

Emerald



Specifications:

Product Code	19055J	1923J	1924J
Seat Width	410 mm	460 mm	510 mm
Width Between Arms	440 mm	465 mm	540 mm
Overall Width	490 mm	520 mm	610 mm
Anterior Seat Height	500 - 590 mm	500 - 590 mm	500 - 590 mm
Posterior Seat Height	475 - 565 mm	475 - 565 mm	475 - 565 mm
SWL	200 kg	200 kg	200 kg