


Cura 1 Premium Bed Pad

Product Code: 92626

Raise an alarm when pressure is removed from the pad. For use with cura1 Premium falls monitors.



Product Features:

- Simple operation and maintenance
- Affordable price
- Easy to clean
- For use with cura1 Premium falls monitors
- USB break cord safety feature
- All bed / chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress
- Warranty: 12 months

Specifications:

Product Code	92626
Length	760 mm
Width	500 mm
SWL	30 - 200 kg