

Alrick 2400 Acute Care Bed

Product Code: 7-24

The 2400 Acute Care Bed is engineered with a robust design to ensure safety and stability. It boasts an impressive safe working load capacity of 350kg to 400kg.



Specifications:

| Size | King Single | Extra Wide | Double | Queen Double |
|---------------------|----------------|----------------|----------------|----------------|
| External Dimensions | 1165 x 2210 mm | 1295 x 2210 mm | 1465 x 2210 mm | 1595 x 2210 mm |
| Sleeping Surface | 1070 x 2000 mm | 1200 x 2000 mm | 1370 x 2000 mm | 1500 x 2000 mm |
| Hi-Lo Position | 390 - 830 mm | 390 - 830 mm | 390 - 830 mm | 390 - 830 mm |
| Underbed Clearance | 160 mm | 160 mm | 160 mm | 160 mm |
| Back / Knee Angle | 77° / 40° | 77° / 40° | 77° / 40° | 77° / 40° |
| TR / ATR Angle | 16° / 16° | 16° / 16° | 16° / 16° | 16° / 16° |
| SWL | 400 kg | 400 kg | 350 kg | 350 kg |

Ensuring safety and stability in an Acute Care environment.



Product Features:

- Australian Made and Owned
- Vertical Lift
- Electric Auto-Regression Backrest
- CPR Backrest
- Electric Kneebreak
- Manual Vascular Support
- 350 - 400kg Safe Working Load
- Back Lit Handset
- Battery Back Up
- Under Bed Light
- Ironstone Powder Coated End Boards with White Laminate Infills
- APO Grey and Ironstone Powder Coated Frame



Forward & Reverse
Trendelenburg



Chair Position
Function



150mm Central
Locking Castors



Bluetooth
Connectivity
***handset must be
ordered separately***



Optional Features - Added at time of manufacture

- Electric Vascular Support
- Pull Out Extension
- Attendant Control
- Split Platform Upgrade
- USB Charging Ports
- Custom Powder Coat Colour
- Custom End Board Colours
- Custom Length of Sleeping Surface (2m+)

Accessories - Added after manufacture

- Mattresses
- Horizontal 400HR Rails (pair)
- Padded Cotside Protectors (pair)
- Self Help Pole
- IV Pole