

# Sentida SC Bed

Product Code: 78082C

Customisable & efficient nursing care bed for demanding environments.



New concepts for living in old age are in demand. These include the desire to be provided with more homeliness, individuality, safety, comfort and flexibility in the equipment of corresponding living arrangements.

With the new nursing bed platform sentida sc and customied concepts, we can respond to the different requirements even more flexibly, efficiently and demand-oriented. We want to support the maintenance and promotion of individual mobility as an overarching and essential goal in care in a solution-oriented manner.

## SafeFree® side rails – safe mobilization

- Optimal for getting out of bed due to its division,
- Individually positionable,
- Easy to grasp and embrace,
- A mobilisation aid for getting up – alone or with support.
- Fear reduction without unnecessary barriers

## Product Features:

- Fully customisable with various models, decor and designs
- SafeFree® side rails – safe mobilization regardless of body size
- MobiStick 2 – manual mobilization assistant
- SafeControl® handset – clear, ergonomic and intuitive
- Double retraction of the lying surface – safe and comfortable mobilisation in bed
- Automatic positions – quick, easy and ergonomic thanks to the one-button philosophy.
- SafeFree® in orientation position – fear reduction without unnecessary barriers
- SafeSense® – autonomous, safe bed exit with individual call notification
- 3-stop strategy – situation- and need-based care, provision and support
- Underbed clearance for mobile patient lifters at any position
- Quality German made construction

## Accessories include:

- Mattresses
- Padded Bed Rail Protectors (pair)
- 200mm Bed Extension
- 200mm Foam Bolster
- IV Pole
- Comfort LED Reading Light
- Lifting Pole with Triangle
- Mobi Stick
- Bed Pole
- SafeLift
- SafeSense



## Specifications:

Product Code	78082C	78256C	78500C
Size	Single	XL	XXL
External Dimensions (max)*	1040 x 2080 mm	1170 x 2080 mm	1370 x 2060 mm
Sleeping Surface	870/900 x 2000 mm	1000 x 2000 mm	1200 x 2000 mm
Back Adjustment (retraction by 120 mm)	70° ± 5°	70° ± 5°	70° ± 5°
Thigh Adjustment (retraction by 50 mm)	30° ± 5°	30° ± 5°	30° ± 5°
Head / Foot Swivel	17° / 17°	17° / 17°	17° / 17°
Hi-Lo Position depending on castor)	230 – 800 mm	230 – 800 mm	230 – 800 mm
Bed Weight (depending on model)	115 kg	130 kg	160 kg
Max Patient Weight	235 kg	235 kg	265 kg
SWL	270 kg	270 kg	300 kg

\*External dimensions differ depending on the design – please contact the team for details.

Issue Date: 030122  
Revision Date: 231225 | V4