

Staebel Standard Mobile Overbed Table

Product Code: 4215

A very popular mobile table with a number of useful features, making it an ideal choice in many situations. With its easy lever operated height adjustment and smooth curved U-design, it fulfils all needs.



One of our most versatile and sturdy over bed tables.

A very popular mobile table with a number of useful features, making it an ideal choice in many situations. With its easy lever operated height adjustment and smooth curved U-design, it fulfils all needs.

Product Features:

- Spring assisted lever height adjustment – lever stays with top eliminating need to bend to adjust table height
- Powdercoated U-shape base, round tube minimises injury risk
- Swivel castors, two are lockable



- 1** **Inverted lift** column to ensure liquids wont leak inside in the event of spills
- 2** Recessed top with wrap to **avoid exposed** seams for infection control
- 3** Spring assist raise and lower making operation **one-handed and extra easy** for staff.
- 4** **Round base** to ensure patients can rest their feet on their without risk of skin tears
- 5** **Additional** locking nut for support if required



Specifications:

Product Code	4215
Top Width	450 mm
Top Length	750 mm
Load Capacity	20 kg

