

Alrick 2300 Tilt Companion Bed

Product Code: 7-23COM

The Alrick 2300 Tilt Companion Bed provides comfort and convenience for users and caregivers alike, designed to meet their needs effectively. This innovative design incorporates two independent bed frames, maximising flexibility and enhancing care outcomes. Users have the option to choose electric operation for both sides or a static frame on one side, depending on care needs.



The beds are connected with a bracket system, allowing them to be used together or separated effortlessly as required.

The electric bed configuration offers adjustable height, backrest, and kneebreak functions, while the static option provides a stable sleeping surface supported by 125mm castors for easy manoeuvrability. Both sides have an independent Safe Working Load of 300kg, ensuring robust support and safety. Available in King Double width, the Alrick 2300 Tilt Companion Bed can be fully customised to meet individual preferences.



Forward & Reverse Trendelenburg



Chair Position



125mm Braking & Directional Castors



Bluetooth Connectivity
handset must be ordered separately

Product Features:

- Australian Made and Owned
- Vertical Lift
- Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- 1 x Electric Bed & 1 x Fixed Height Bed
- Beds can be easily separated
- 300kg SWL each side
- Pull Out Extension (both sides of the beds)
- Back Lit Handset
- Under Bed Light
- Curved Natural Oak End Boards
- APO Grey and Ironstone Powdercoated Frame
- Forward and Reverse Trendelenburg
- Chair Position
- 125mm Braking & Directional Castors
- Bluetooth Connectivity
*handset must be purchased separately *

High flexibility!

As the bed is securely connected by a **joining bar in between the two beds**, it can be conveniently split into two separate beds. This feature not only facilitates effortless delivery and transportation into and out of rooms, but also allows for future versatility in configuration and use.

Optional Features – Added at time of manufacture

- 2 x Electric Beds
- CPR Backrest
- Electric Vascular Support
- Central Locking Castors
- Bluetooth Handset (HC40)
- Battery Back Up
- USB Charging Ports
- Under Bed Light with Custom Colour
- Custom Powdercoat Colours
- Custom End Board Styles and Colours
- Custom Length Sleeping Surface (1.5m+)
- Padded York Surround

Accessories – Added after manufacture

- Mattresses
- Horizontal 400HR Rails (pair)
- Padded Cotside Protectors (pair)
- Clamp on Support Rail
- StandEzy Straight Handle
- StandEzy Hook Handle
- Self Help Pole
- 200mm Foam Bolster
- Bluetooth Handset (7-BT-4-HB200)
- Emergency Bed Power Unit
- IV Pole
- Combo Grab Handle
- T-Bar Buffer
- Reading Lamp

Specifications:

Type	Two Electrics	One Electric, One Fixed
Size	King Double	King Double
External Dimensions (Overall)	1840 x 2050 mm	1840 x 2050 mm
External Dimensions (Each)	920 x 2050 mm	920 x 2050 mm
Sleeping Surface	900 x 2000 mm	900 x 2000 mm
Hi-Lo Position	360 – 830 mm	400 mm
Underbed Clearance	170 mm	140 mm
Backrest / Knee Angle	77° / 40°	-
TR / ATR Angle	16° / 16°	-
SWL	300 kg (per side)	300 kg (per side)

*Standard specifications listed above. For more details regarding width sizes, refer to the script form.

Issue Date: 240424
Revision Date: 280825 | V2

