

Alrick 2300 Tilt Care Bed

Product Code: 7-23T

The Alrick 2300 Tilt Care Bed is the most versatile Hi-Low and tilt bed available for the healthcare industry. Its robust frame supports up to 300kg SWL, accommodating users with diverse weight requirements. Available in sizes ranging from single to king double width, this bed is equally suitable for both in-home and facility care settings, combining aesthetic appeal with practicality for users and carers alike.



Product Features:

- Australian Made and Owned
- Vertical Lift
- Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- 300kg SWL
- Pull Out Extension
- Back Lit Handset
- Under Bed Light (Warm White)
- Curved Natural Oak End Boards
- APO Grey and Ironstone Powder Coated Frame

Optional Features - Added at time of manufacture

- Partner Bed Set Up
- CPR Backrest
- Electric Vascular Support
- Central Locking Castors
- Bluetooth Handset (HC40)
- Attendant Control
- USB Charging Ports
- Underbed Light with Custom Colour
- Battery Back Up
- Custom Powdercoat Colour
- Custom Endboard Style and Colour
- Custom Length Sleeping Surface (1.5m+)

Accessories - Added after manufacture

- Mattresses
- Horizontal Rails
- Padded Cotside Protector
- Clamp On Support Rail
- StandEzy Straight Handle
- StandEzy Hook Handle
- Self Help Poles
- 200mm Foam Bolster
- Bluetooth Handset (7-BT-4-HB200)
- Emergency Bed Power Unit
- IV Pole
- Combo Grab Handle
- T-Bar Buffer
- Reading Lamp



Forward & Reverse Trendelenburg



Chair Position



125mm Braking & Directional Castors



Bluetooth Connectivity
handset must be ordered separately

The most versatile tilt bed for all healthcare environments.




Pull Out extension
as standard

Stocked Endboard Styles & Colours

Custom Endboard styles & colours available upon request.



Laminate
Curved Style



Upholstered
Jersey Headboard &
Plain Footboard



Timber
Eltham Style



Natural Oak
Laminate



Ocean
Upholstered




Ash
Upholstered



Burnt Walnut
Timber




Beech
Timber



Saddle
Timber



Custom
Endboards



Specifications:

Size	Single	King Single	Extra Wide	Double	Queen Double	King Double
External Dimensions	920 x 2040 mm	1090 x 2040 mm	1220 x 2040 mm	1390 x 2040 mm	1520 x 2040 mm	1820 x 2040 mm
Sleeping Surface	900 x 2000 mm	1070 x 2000 mm	1200 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm	1800 x 2000 mm
Hi - Lo Position	360 - 800 mm	360 - 800 mm	360 - 800 mm	360 - 800 mm	360 - 800 mm	360 - 800 mm
Underbed Clearance	140 mm	140 mm	140 mm	140 mm	140 mm	140 mm
Back / Knee Angle	77° / 40°	77° / 40°	77° / 40°	77° / 40°	77° / 40°	77° / 40°
TR / ATR Angle	16° / 16°	16° / 16°	16° / 16°	16° / 16°	16° / 16°	16° / 16°
SWL	300 kg	300 kg	300 kg	300 kg	300 kg	300 kg

Issue Date: 050624
Revision Date: 280825 | V2