

Alrick 2300 Partner Bed

Product Code: 7-23TP

Specifically designed to cater to the individual needs of two users while maintaining a connected and cohesive design through integrated endboards. This bed features electric hi-low functionality that operates in unison, ensuring synchronised adjustments across both sides. The bed includes individual electric backrests and kneebreaks, each capable of supporting a safe working load of 200kg.



Product Features:

- Australian Made and Owned
- Vertical Lift
- Auto-Regression Backrests
- Electric Kneebreak
- Manual Vascular Support
- Easily Separated
- 200kg SWL (each side)
- 2 x Electric Beds joined by End Boards
- Pull Out Extensions
- Back Lit Handsets
- Underbed Lights (Warm White)
- Curved Natural Oak End Boards
- APO Grey and Ironstone Powder Coated Frame



Forward & Reverse
Trendelenburg



Chair Position



125mm Braking
& Directional
Castors



Bluetooth
Connectivity
***handset must be
ordered separately***

Specifications:

Size	King Double
External Dimensions	1830 x 2050 mm
Sleeping Surface	1800 x 2000 mm
Hi - Lo Position	360 - 800 mm
Underbed Clearance	140 mm
Back / Knee Angle	77° / 40°
TR / ATR Angle	16° / 16°
SWL	200 kg (per side)

Issue Date: 260424
Revision Date: 2300326 | V3



The bed is connected by shared endboards, which allows the bed to be split into two separate beds for versatility and transport.

Optional Features – Added at time of manufacture

- Floorline Available
- CPR Backrest
- Electric Vascular Support
- Central Locking Castors
- Bluetooth Handset (HC40)
- Battery Back Up
- USB Charging Ports
- Underbed Light with Custom Colours
- Custom Powder Coat Colours
- Custom Endboard Styles and Colours
- Custom Length Sleeping Surface (1.5m+)
- Padded York Surround



Pull Out extension as standard

Stocked Endboard Styles & Colours


Custom Endboard styles & colours available upon request.



Laminate
Curved Style



Upholstered
Jersey Headboard &
Plain Footboard



Timber
Eltham Style



Natural Oak
Laminate



Ocean
Upholstered




Ash
Upholstered



Burnt Walnut
Timber



Beech
Timber



Saddle
Timber



Accessories – Added after manufacture

- Mattresses
- Horizontal Rails (pairs)
- Padded Cotside Protectors
- Clamp on Support Rail
- StandEzy Straight Handle
- StandEzy Hook Handle
- Self Help Poles
- 200mm Foam Bolster
- Bluetooth Handset (7-BT-4-HB200-2)
- UPS Emergency Bed Power Unit
- IV Pole
- Combo Grab Handle
- T-Bar Buffer
- Reading Lamp

