

Alrick 2300 Split Bed

Product Code: 7-23TS

The Alrick 2300 Split Bed is designed to cater separately to the needs of two users within a single bed frame, offering a split sleeping surface for individual comfort. Built with heavy-duty actuators, this innovative bed includes combined electric hi-low functionality, along with individual electric adjustments for backrests and kneebreaks, all supporting a total safe working load of 300kg.



Product Features:


- Australian Made and Owned
- Vertical Lift
- 2 x Independent Electric Backrest
- 2 x Independent Electric Kneebreak
- Chair Position
- Manual Vascular Support
- Under Bed Light (Warm White)
- Back Lit Handset
- Battery Back Up
- Curved Natural Oak End Boards
- APO Grey and Ironstone Powder Coated Frame



Forward & Reverse
Trendelenburg



**Pull Out
Extension**



**125mm Twin
Breaking Castors**



**Bluetooth
Connectivity**
***handset must be
ordered separately***

Specifications:

Size	Double	Queen Double
External Dimensions	1390 x 2050 mm	1520 x 2050 mm
Sleeping Surface	1370 x 2000 mm	1500 x 2000 mm
Low Position	360 mm	390 mm
High Position	800 mm	830 mm
Underbed Clearance	150 mm	170 mm
Back / Knee Angle	77° / 40°	77° / 40°
TR / ATR Angle	16° / 16°	16° / 16°
SWL	300 kg	300 kg



The most versatile split bed for all healthcare environments.

Optional Features – Added at time of manufacture

- CPR Backrest
- Electric Vascular Support
- Central Locking Castors
- Bluetooth Handset
- Underbed Light with Custom Colour
- USB Charging Ports
- Custom Powder Coat Colours
- Custom End Board Styles and Colour
- Custom Length Sleeping Surface (1.5m+)



Stocked Endboard Styles & Colours

Custom Endboard styles & colours available upon request.



Laminate
Curved Style



Upholstered
Jersey Headboard &
Plain Footboard



Timber
Eltham Style



Natural Oak
Laminate



Ocean
Upholstered




Ash
Upholstered



Burnt Walnut
Timber



Beech
Timber




Saddle
Timber



Accessories – Added after manufacture

- Mattresses
- Horizontal 400HR Rails (pair)
- Padded Cotside Protectors (pair)
- Clamp on Support Rail
- StandEzy Straight Handle
- StandEzy Hook Handle
- Self Help Poles
- Bluetooth Handset (7-BT-4-HB200-2)
- 200mm Foam Bolster
- Emergency Bed Power Unit
- IV Pole
- Combo Grab Handle
- T-Bar Buffer
- Reading Lamp



Issue Date: 260424
Revision Date: 2300326 | V3