

Alrick Overbed Table

Product Code: 7-A6010

The Alrick Overbed Table fully compliant with Australian standards. It has a height adjustable design with a C base shape to allow for easy access. It also has 4 Castors for easy manoeuvrability.



Begin with a standard style for our Overbed Tables, or branch out and talk with us about a customised product. We use moisture-resistant laminate, and timber tops are coated in a bio-guard coating. All our products in the range are fully compliant with Australian standards. Our products are chlorine-free and environmentally friendly, and they're available in various colours and designs.



Product Features:

- Australian Made & Owned
- Easy to adjust the height
- Powdercoated C-shape base, round tube minimises injury risk
- Swivel castors, two are lockable
- Available in Flat Laminated Top or Flat Timber Top
- Standard Laminate or Timber stain options available, with the option to choose a custom colour



Specifications:

Product Code	7-A6010
Height	755 - 1080 mm
Top Width	390 mm
Top Length	790 mm
Load Capacity	20 kg

Timber Colour Options:



Laminate Colour Options :

