

Horizon Acute Floorline Bed

Product Code: 7-AEN7

The Horizon Acute Floorline Bed is one of the lowest beds available on the market and suitable for acute care geriatric and falls risk environments. The bed lowers to floor level to prevent injuries or accidents caused by falls.




Product Features:

- Floorline Model
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- SWL 200kg
- Back Lit Handset
- CPR Release
- Battery Back Up
- Plastic Endboards
- Ironstone Powder Coated Frame as standard



Forward & Reverse
Trendelenburg



Chair
Position



100mm Twin
Wheel Braking
Castors



Bluetooth
Connectivity
***handset must be
ordered separately***

Specifications:

Size	Single
External Dimensions	920 x 2115 mm
Sleeping Surface	900 x 2000 mm
Low Position	135 mm
High Position	750 mm
Backrest / Knee Angle	62° / 35°
TR / ATR Angle	15° / 15°
SWL	200 kg



The universal floorline bed that suits low care, high care & bariatric care.

Optional Features – Added at time of manufacture

- Underbed Light
- Bluetooth Handset (HC40)
- USB Charging Ports



White Plastic Endboards

Accessories – Added after manufacture

- 400HR Horizontal Rails
- 300kg Handset Upgrade
- 200mm Bed Extension
- Mattresses
- StandEzy Hook Handle
- Self Help Pole
- UPS Emergency Power Bed Unit
- IV Pole
- Combo Grab Handle
- EN T-Bar Buffer
- Reading Lamp

