

# Alrick Cira Tilt Bed

Product Code: 7-CIRA

Offers refined functionality with a homely feel. Designed for more complex needs, this bed supports bariatric or high-dependency care without compromising on comfort or aesthetics. It's a reassuring presence in any care setting where dignity and durability are paramount.



## Specifications:

| Size                | Single         | King Single    | Extra Wide     | Double         | Queen Double   | King Double    |
|---------------------|----------------|----------------|----------------|----------------|----------------|----------------|
| External Dimensions | 1050 x 2060 mm | 1220 x 2060 mm | 1350 x 2060 mm | 1520 x 2060 mm | 1650 x 2060 mm | 1950 x 2060 mm |
| Sleeping Surface    | 900 x 2000 mm  | 1070 x 2000 mm | 1200 x 2000 mm | 1370 x 2000 mm | 1500 x 2000 mm | 1800 x 2000 mm |
| Hi - Lo Position    | 360 - 800 mm   | 360 - 800 mm   | 360 - 800 mm   | 360 - 800 mm   | 360 - 800 mm   | 360 - 800 mm   |
| Underbed Clearance  | 140 mm         | 140 mm         | 140 mm         | 140 mm         | 140 mm         | 140 mm         |
| Back / Knee Angle   | 77° / 40°      | 77° / 40°      | 77° / 40°      | 77° / 40°      | 77° / 40°      | 77° / 40°      |
| TR / ATR Angle      | 16° / 16°      | 16° / 16°      | 16° / 16°      | 16° / 16°      | 16° / 16°      | 16° / 16°      |
| SWL                 | 300 kg         | 300 kg         | 300 kg         | 300 kg         | 300 kg         | 300 kg         |



## Product Features:

- Australian Made and Owned
- Four-part mobilisation rails
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- SWL 300kg
- Manual Vascular Support
- Back Lit Handset
- Underbed Light (Warm White)
- 25mm MDF Laminate Endboards
- APO Grey and Ironstone Powder Coated Frame



Forward & Reverse  
Trendelenburg



Chair  
Position



125mm Braking  
& Directional  
Castors



Bluetooth  
Connectivity  
\*\*\*handset must be  
ordered separately\*\*\*

## Stocked Endboard Colours & Styles

*Custom Endboard, styles and colours also available.*



### Laminate

Square Headboard & Footboard



**Natural Oak**  
Laminate



**Ocean**  
Upholstered



**Ash**  
Upholstered



### Upholstered

Plain Headboard & Footboard



## Optional Features:

- Central Locking Castors
- Bluetooth Handset (HC40)
- Attendant Control
- Battery Back Up
- USB Charging Ports
- Custom Underbed Light Colour
- Custom Powdercoat Colour
- Custom Endboard Style and Colour

## Accessories - Added after manufacture

- Bluetooth Handset (7-BT-4-HB200)
- Mattresses
- Self Help Pole
- Emergency Power Bed Unit
- IV Pole
- T-Bar Buffer
- Reading Lamp
- Bed Rail Protectors (3983)

