

Staebel EN7 Floorline Bed

Product Code: 7-EN7

The Staebel EN7 Floorline Bed is one of the lowest beds available on the market and designed for intensive use. The bed lowers to floor level to prevent injuries and accidents caused by falls from beds. All actuators and mechanisms are concealed beneath the platform, allowing the bed to be constructed with no end pedestals!




Product Features:

- Floorline Model
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- SWL 200kg
- Back Lit Handset
- Natural Oak Slender Arch End Boards
- Ironstone Powder Coated Frame as standard



Forward & Reverse
Trendelenburg



Chair Position
Function



75mm Twin
Wheel Braking
Castors



Bluetooth
Connectivity
***handset must be
ordered separately***

Specifications:

Size	Single	Wide	King Single
External Dimensions	920 x 2115 mm	1020 x 2115 mm	1090 x 2115 mm
Sleeping Surface	900 x 2000 mm	1000 x 2000 mm	1070 x 2000 mm
Hi-Lo Position	100 - 725 mm	100 - 725 mm	100 - 725 mm
Backrest / Knee Angle	62° / 35°	62° / 35°	62° / 35°
TR / ATR Angle	15° / 15°	15° / 15°	15° / 15°
SWL	200 kg	200 kg	200 kg

Issue Date: 281022
Revision Date: 300626 | V4


The universal floorline bed that suits low care, high care & bariatric care.

Optional Features – Added at time of manufacture


- Underbed Light (Warm White)
- Underbed Light with Custom Colours
- Bluetooth Handset (HC40)
- Battery Back Up
- USB Charging Ports
- Dual Locking Castors
- 100mm Twin Wheel Castors
- Attendant Control
- CPR Backrest
- Custom End Board Styles and Colours

Stocked Endboard Styles & Colours

Custom Endboard styles & colours available upon request.



Laminate
Slender Arch Style



Upholstered
Jersey Headboard &
Plain Footboard




Timber
Eltham Style



Natural Oak
Laminate



Ocean
Upholstered



Ash
Upholstered



Burnt Walnut
Timber



Beech
Timber



Saddle
Timber

Accessories – Added after manufacture

- 300kg Handset (Single & King Single Only)
- 200mm Bed Extension
- Mattresses
- Horizontal Bed Rails
- StandEzy Hook Handle
- Self Help Pole
- IV Pole
- Combo Grab Handle
- EN T-Bar Buffer
- Reading Lamp
- Bluetooth Handset (7-BT-4-HB200)
- UPS Emergency Power Bed Unit

