

Horizon Kova Acute Bed

Product Code: 7-KOVA

The Kova Acute Bed is engineered with a robust, durable design to ensure safety and stability. It boasts an impressive safe working load capacity of 500kg.



Specifications:

Size	King Single	Extra Wide	Double	Queen Double
External Dimensions	1165 x 2210 mm	1295 x 2210 mm	1465 x 2210 mm	1595 x 2210 mm
Sleeping Surface	1070 x 2000 mm	1200 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm
Hi-Lo Position	390 - 740 mm	390 - 740 mm	390 - 740 mm	390 - 740 mm
Underbed Clearance	160 mm	160 mm	160 mm	160 mm
Back / Knee Angle	68° / 40°	68° / 40°	68° / 40°	68° / 40°
TR / ATR Angle	16° / 16°	16° / 16°	16° / 16°	16° / 16°
Max Patient Weight	400 kg	400 kg	400 kg	400 kg
SWL	500 kg	500 kg	500 kg	500 kg

Ensuring the highest safety and stability in all healthcare environments.



Product Features:

- Australian Made and Owned
- Vertical Lift
- Electric Auto-Regression Backrest
- CPR Backrest
- Electric Kneebreak
- Manual Vascular Support
- 500kg Safe Working Load
- Back Lit Handset
- Battery Back Up
- Under Bed Light (Warm White)
- Ironstone Powder Coated End Boards with White Laminate Infills
- APO Grey and Ironstone Powder Coated Frame



Forward & Reverse
Trendelenburg



**Chair
Position**



**150mm Central
Locking Castors**



**Bluetooth
Connectivity**
***handset must be
ordered separately***



Optional Features – Added at time of manufacture

- Electric Vascular Support
- Pull Out Extension
- Attendant Control
- Split Platform Upgrade
- USB Charging Ports
- Under Bed Light with Custom Colours
- Custom Powder Coat Colour
- Custom End Board Colours
- Custom Length of Sleeping Surface (2m+)

Accessories – Added after manufacture

- Mattresses
- Horizontal 400HR Rails (pair)
- Padded Cotside Protectors (pair)
- Self Help Pole
- IV Pole