

# Alrick Lumi Bed




Product Code: 7-LUMI

Offering refined functionality with a homely feel. With smooth adjustability, safety-focused features, and a warm, non-institutional aesthetic, the Lumi is ideal for aged care homes and assisted living environments that want to blend comfort with care.



## Specifications:

| Size                  | King Single    | Double         | Queen Double   |
|-----------------------|----------------|----------------|----------------|
| External Dimensions   | 1220 x 2160 mm | 1520 x 2160 mm | 1650 x 2160 mm |
| Sleeping Surface      | 1070 x 2000 mm | 1370 x 2000 mm | 1500 x 2000 mm |
| Hi-Lo Position        | 200 - 725 mm   | 200 - 725 mm   | 200 - 725 mm   |
| Backrest / Knee Angle | 62° / 35°      | 62° / 35°      | 62° / 35°      |
| TR / ATR Angle        | 15° / 15°      | 15° / 15°      | 15° / 15°      |
| SWL                   | 250 kg         | 250 kg         | 250 kg         |




## Product Features:

- Australian Made and Owned
- Four-part mobilisation rails
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- SWL 250kg
- Back Lit Handset
- Natural Oak Laminate Square Endboards
- Ironstone Powdercoated Frame as standard



Forward & Reverse  
**Trendelenburg**



**Chair Position**



**75mm Twin**  
Wheel Braking  
Castors



**Bluetooth Connectivity**  
\*\*\*handset must be ordered separately\*\*\*




## Stocked Endboard Colours & Styles

*Custom Endboard, styles and colours also available.*



**Laminated**  
Square Headboard & Footboard



**Upholstered**  
Plain Headboard & Footboard



**Natural Oak**  
Laminate



**Ocean**  
Upholstered



**Ash**  
Upholstered



## Optional Features:

- Dual Locking Castors
- Bluetooth Handset (HC40)
- Attendant Control
- Battery Back Up
- Two or Four Rails Mobilisation rails
- Underbed Lighting (Warm White)
- Underbed Lighting with Custom Colours
- USB Charging Ports
- Custom Endboard Style and Colour



## Accessories - Added after manufacture

- Bluetooth Handset (7-BT-4-HB200)
- Mattresses
- Self Help Pole
- Emergency Power Bed Unit
- IV Pole
- EN T-Bar Buffer
- Reading Lamp
- Bed Rail Protector (3983)



Issue Date: 280825  
Revision Date: 090726 | V3