

Alrick Lumi Bed


Product Code: 7-LUMI

Offering refined functionality with a homely feel. With smooth adjustability, safety-focused features, and a warm, non-institutional aesthetic, the Lumi is ideal for aged care homes and assisted living environments that want to blend comfort with care.



Specifications:

Size	King Single	Double	Queen Double
External Dimensions	1220 x 2160 mm	1520 x 2160 mm	1650 x 2160 mm
Sleeping Surface	1070 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm
Hi-Lo Position	200 - 725 mm	200 - 725 mm	200 - 725 mm
Backrest / Knee Angle	62° / 35°	62° / 35°	62° / 35°
TR / ATR Angle	15° / 15°	15° / 15°	15° / 15°
SWL	250 kg	250 kg	250 kg




Product Features:

- Australian Made and Owned
- Four-part mobilisation rails
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- SWL 250kg
- Back Lit Handset
- 25mm MDF Laminate Endboards
- Ironstone Powdercoated Frame as standard



Forward & Reverse Trendelenburg



Chair Position



75mm Twin Wheel Braking Castors



Bluetooth Connectivity
handset must be ordered separately




Stocked Endboard Colours & Styles

Custom Endboard, styles and colours also available.



Laminate
Square Headboard & Footboard



Upholstered
Plain Headboard & Footboard



Natural Oak
Laminate



Ocean
Upholstered



Ash
Upholstered



Optional Features:

- Dual Locking Castors
- Bluetooth Handset (HC40)
- Attendant Control
- Battery Back Up
- Underbed Lighting (Warm White)
- Underbed Lighting with Custom Colours
- USB Charging Ports
- Custom Endboard Style and Colour



Accessories - Added after manufacture

- Bluetooth Handset (7-BT-4-HB200)
- Mattresses
- Self Help Pole
- Emergency Power Bed Unit
- IV Pole
- EN T-Bar Buffer
- Reading Lamp
- Bed Rail Protector (3983)

