

Alrick Luxa Homecare Bed

Product Code: 7-LUXA

A premium homecare bed that is designed to give the user ultimate comfort and care, while maintaining dignity and style in the bedroom. This bed features electric auto-regression backrest, manual vascular support, a high height of 725mm and a low bed height of 150mm.



Product Features:

- Australian Made and Owned
- Floorline Bed
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- 250kg SWL
- Back Lit Handset
- Battery Back Up
- Ironstone Powder Coated Frame
- Full Surround Design



Forward & Reverse
Trendelenburg



Chair
Position





75mm Twin
Wheel Braking
Castors



Bluetooth
Connectivity
***handset must be
ordered separately***

Optional Features – Added at time of manufacture

- Dual Locking Castors
- Bluetooth Handset (HC40)
- Under Bed Light (Warm White)
- Under Bed Light with Custom Colours
- USB Charging Ports
- Timber Surround
- Custom Powder Coat Colours
- Custom End Board Styles and Colours
- Custom Sleeping Surface Length (2m+)
- Partner Bed Set Up



Stocked Styles & Colours

Custom styles & colours available upon request.



Upholstered

Plain Padded Headboard & Plain Full Surround



Upholstered

Jersey Padded Headboard & Plain Full Surround



Foot section of surround available in Low or High Profile.



Natural Oak
Laminate



Ocean
Upholstered



Ash
Upholstered

Accessories - Added after manufacture

- Mattresses
- Bluetooth Handset (7-BT-4-HB200)
- Luxa Horizontal 400HR Bed Rails (pair)
- Padded Cotside Protector (pair)
- Luxa Clamp on Support Rail
- Luxa 200mm Bed Extension
- 200mm Foam Bolster (order with Bed Extension)
- Emergency Bed Power Unit
- IV Pole
- EN T-Bar buffer
- Reading Lamp



Specifications:

Size	Single	King Single	Double	Queen Double
Overall Dimensions	1000 x 2160 mm	1170 x 2160 mm	1470 x 2160 mm	1600 x 2160 mm
Sleeping Surface	900 x 2030 mm	1070 x 2030 mm	1370 x 2030 mm	1500 x 2030 mm
Hi - Lo Position	150 - 725 mm	150 - 725 mm	170 - 725 mm	170 - 725 mm
TR / ATR Angle	62° / 35°	62° / 35°	62° / 35°	62° / 35°
Back / Knee Angle	15° / 15°	15° / 15°	15° / 15°	15° / 15°
SWL	250 kg	250 kg	250 kg	250 kg

Issue Date: 130224
Revision Date: 230326 | V3