

Nestie Novus Maternity Partner Bed

Product Code: 7-NOVUS

The Nestie Novus Maternity Partner Bed allows parents to stay together and care for their newborn, providing a safe space and providing ease for care staff.



Specifications:

Size	Double	Queen Double
Overall Dimensions	1465 x 2040 mm	1595 x 2040 mm
Sleeping Surface	1370 x 2000 mm	1500 x 2000 mm
Hi-Lo Range	390 – 830 mm	390 – 830 mm
Underbed Clearance	160 mm	160 mm
Back / Knee Angle	77° / 40°	77° / 40°
TR / ATR Angle	16° / 16°	16° / 16°
SWL	300 kg	300 kg

Issue Date: 160425
Revision Date: 270426 | V3

The Nestie Novus Maternity Partner Bed enables parents to stay together and remain actively involved in caring for their new born baby in a nurturing environment, enhancing the overall childbirth experience. With height adjustability empowering nurses to provide exceptional care and minimising strain. The removable leg stirrups which provides a safe, comfortable space, enabling post birth suturing and medical check-ups with ease.

Product Features:

- Australian Made and Owned
- Electric Hi-Low
- Electric Backrest
- Electric Kneebreak
- Battery Backup
- Underbed Light (Warm White)
- Backlit Handset
- Safe Working Load 300kg
- Ironstone steel frame and Matt Classic White Infill Endboards
- APO Grey and Ironstone Powdercoated Frame
- Birthing Stirrup Bar *(Stirrups to be purchased separately)*



Forward & Reverse
Trendelenburg



Manual CPR
Release



**150mm Twin
Wheel** Central
Locking Castors



**Bluetooth
Connectivity**
***handset must be
ordered separately***



Optional Features – Added at time of manufacture

- Upgrade to 400kg SWL
- Underbed Light with Custom Colours
- Electric Vascular Support
- Attendant Control
- Bluetooth Handset
- Laminate or Upholstered Endboards

Accessories – Added after manufacture

- Stirrups
- Mattresses
- IV Pole
- Self Help Pole & Triangle
- Bariatric Self Help Pole & Triangle
- 400HR Horizontal Deluxe Rails



*Custom Endboards Shown in
Denim Panache*